



How to Lock the Buttons

We strongly suggest that you keep buttons locked when you have completed setting the time, date, and alarms to prevent accidental re-programming.

- 1) Press **MODE** repeatedly to advance to Time Mode.
- 2) Press and hold **MODE** and **RESET** together at the same time to lock the buttons. A "LOCK" icon will appear in the top left-hand corner of the display. When buttons are locked, only the **LIGHT** button will be working.
- 3) Press and hold **MODE** and **RESET** together at the same time to remove the lock-out feature ("LOCK" icon will disappear).

How to Erase an Alarm

- 1) Press **MODE** repeatedly to advance to Alarm Mode. "AL" will appear at the top of the display.
- 2) Press **RESET** repeatedly to get to specific alarm you want to erase. The alarm number is shown in the upper right-hand corner.
- 3) Press and hold **RESET** until , , etc. disappears from the middle of the display.
- 4) Press **MODE** until you return to Time Mode. The day of the week, current time, and date should now appear. Your alarm(s) have been erased.

How to Change the Battery

The e-pill CADEX VibraPlus runs on a **CR2032 button cell**. At the first sign of low power (dim display), have the battery replaced. We strongly recommend that you take your e-pill Watch to a jeweler for battery replacement, instead of changing it on your own.

CAUTION: Dead battery, replacement of battery, malfunction or repairs will cause all data stored in the memory of the watch to be deleted. It is recommended that you keep separate written copies of all important data. In no event will e-pill, LLC be liable for loss of data or damages arising of such loss.

IMPORTANT INFORMATION

Water Resistance

This watch is **NOT WATER PROOF**. It will withstand splashes and rain but not swimming, bathing, and showering.

DO NOT USE IN SHOWER OR SWIMMING POOL

Care of Your Watch

Avoid exposing the watch to temperature extremes. Avoid dropping the watch. Avoid exposing the watch to gasoline, cleaning solvents, adhesive agents, paint, or aerosol sprays.
DO NOT IMMERSE IN WATER

e-pill® CADEX® VibraPlus™


8 Alarm Vibrating Medication Reminder Watch

Thank you for purchasing the e-pill CADEX VibraPlus 8 Alarm Watch. Your e-pill Watch has the following functions:

- Up to 8 Daily Auto-Repeating Alarms
- Daily Alarms Sound, Vibrate, or do Both
(Alarms are 15 seconds long)
- Lock Feature to Prevent Accidental Re-Programming
- Easy-to-Use Interface with Backlight
- Countdown Timer with Time Limit
(Ex. Alarm from 8:00 AM to 9:00 PM)

LIGHT

ST/STP



MODE

RESET

To change the time from 24 Hour Format to AM/PM Format, Press ST/STP when in the Time Mode.

WARRANTY

The inner components of this watch are warranted to be free from defects in materials and manufacturer's workmanship for one (1) year from the original date of purchase. This warranty does not cover the band, buckle, battery, case, buttons, or crystal. Do not wear this watch in the shower or swimming pool, it is **NOT WATER PROOF**. This warranty expressly excludes incidental and consequential damage caused by use, misuse, neglect, accident, unauthorized modification, or any defects resulting from service by anyone other than the manufacturer. During the warranty period, and upon proof of purchase, this watch will be repaired or replaced (with the same or a similar model) at the manufacturer's option.

Contact: e-pill Medication Reminders (or your distributor)



Web site: www.TabTimer.com.au /atch.com
Phone: 1-88 www.TabTimer.com.au (ernational)


1300 TAB TIMER (1300 822 846)

ABS. 99 127 415 246

The e-pill CADEX VibraPlus watch is only a supplement to any memory system that you currently use to remind yourself to take your medications or perform medical tasks. The user maintains full responsibility to administer medications or perform medical tasks at the appropriate times.

e-pill, LLC assumes no responsibility for any incidental or consequential damages, including but without limitation, damages or injuries resulting from malfunction of the product, mathematical inaccuracy or loss of stored data.


Four Display Modes



Flip through the display modes by pressing the MODE button repeatedly.

How to Set the Time and Date

- 1) Select the Time Mode by pressing **MODE** repeatedly if not already shown. Press and hold **LIGHT** for five (5) seconds. This will cause the seconds to begin flashing.
- 2) Press **MODE** and the minutes will flash. Press **ST/STP** to advance the minutes. Press **MODE** to record the minutes. The hour will now flash.
- 3) Please note the "AM" or "PM" Indicator to the left of the hour digits. Press **ST/SP** to advance the hour, then press **MODE** to record the hour.
- 4) The date will flash. Press **ST/SP** to advance the date, then press **MODE** to record the date.
- 5) The month will flash. Press **ST/SP** to advance the month, then press **MODE** to record the month.
- 6) The year will flash. Press **ST/SP** to advance the year, then press **MODE** to record the year.
- 7) Press **LIGHT** to confirm the Time and Date Settings and return to Time Mode.



To temporarily display the Date in Time Mode, Press **RESET**.

How to Set the Alarm to Beep and/or Vibration

This section explains how to set Beep/Vibration/Both for ALL Alarms, including the Countdown Timer.

- 1) Select the Alarm Mode (AL) by pressing **MODE** repeatedly if not already shown.
- 2) Press **ST/STP** once to turn on the beeping alarm, twice to turn on the vibration, or three times to turn on both alarm and vibration.
 - When sound alarm is on, (())) will display.
 - When vibration is on, \$\$\$ will display.
 - When both sound and vibration are on, (())) \$\$\$ will display.

These symbols can be found above the AM/PM indicator on the left side of the display.



How to Set up to 8 Daily Alarms

- 1) Select the Alarm Mode (AL) by pressing **MODE** repeatedly if not already shown. Alarms are shown in the top right-hand corner, numbered A1-A8.
- 2) To set Alarm 1, press and hold **LIGHT** until hour flashes. Press **ST/STP** to advance the hour, then press **MODE** to record the hour.
- 3) The minutes will flash. Press **ST/STP** to advance the minutes, then press **MODE** to record the minutes.
- 4) Press **LIGHT** to confirm the time entered. A, A, etc. will appear above the time to show that the alarm has been set.
- 5) Press **RESET** to move between alarms (A1-A8). Repeat steps 1-4 for each alarm you want to set.
- 6) When you are finished setting alarms press **MODE** repeatedly until you reach Time Mode.

Alarms auto-reset at midnight and repeat automatically again at the same time the next day.

How to Set the Stopwatch

- 1) Select the Stopwatch Mode (ST) by pressing **MODE** repeatedly if not already shown. "ST" will appear on the top of the display.
- 2) Press **ST/STP** to start/stop the stopwatch.
- 3) Press **RESET** to return to 0.

How to Set the Countdown Interval Timer

The Countdown Timer works within a set time-frame every day - for example, 8:00AM to 9:00PM.

Within these hours the watch will alarm at set intervals - for example, every hour, or every 2 hours. Set in hours and minutes (Minimum Interval 15 min.)

To set the Time Frame/ Time Interval:

- 1) Select the Countdown Timer Mode (TR) by pressing **MODE** repeatedly if it is not already shown.
- 2) Press and hold **LIGHT** until the starting hour flashes. This is when the countdown timer will start every day. Press **ST/STP** to advance to this hour.
- 3) Press **MODE** and the ending hour will flash. This is when the countdown timer will stop every day. This hour must be later than the starting hour. Press **ST/STP** to advance to this hour.

Next Set the Time Interval (Ex: 30 min, 1 hr, 3hrs, etc.):

- 4) Press **MODE** and the hour flash. Press **ST/STP** to advance the hour, then press **MODE** to record the hour.
- 5) The minutes will flash. Press **ST/STP** to advance the minutes, then press **MODE** to record the minutes.
- 6) Press **LIGHT** to confirm the timer settings.

ST/STP turns the Countdown Interval Timer on and off. The timer will automatically restart daily.

Make sure to set the Alarm to Beep and/or Vibration.