## English

**Sonic Alert**  
Specialized Hearing Solutions  
**Model SBT600ss**  
**Sonic Boom Alarm Clock & Bed Shaker**  
with 5-Level Dimmer & USB Charging

### OPERATION INSTRUCTIONS

#### IMPORTANT:

- Please read these instructions carefully before use and retain for future reference. See operating instructions on back.
- Use only the supplied or recommended Sonic Alert adaptor to connect the unit.

#### Warning:




To prevent fire or shock hazard do not expose this appliance to rain or moisture. Before connecting to mains/ power make sure that the requirements specified on the apparatus corresponds to the available power supply in that country i.e. AC100-240V~, 50/60 Hz.

- The ventilation should not be impeded by covering the ventilation openings with items such as, newspaper, table cloth, curtains, etc;
- No naked flame sources, such as lighted candles should be placed on the apparatus.
- It shall not be exposed to dripping or splashing and that no objects filled with liquids, such as vases, shall be placed on the apparatus.
- Do not place the product in closed bookcases or racks without proper ventilation.
- Where the mains/ power plug is used as the disconnect device, the disconnect device shall remain readily operable.

#### Test Mode:

- 1) Press and hold the Snooze and Time buttons for 3 seconds to enter the Test Mode. The display will show "7ES7" to indicate it is in TEST Mode.
- 2) If within 1 minute no button is pressed, then it will exit the Test Mode and return back to the Time Mode.
- 3) Press any button within 1 minute (except the Snooze & Time buttons, which has no effect), it will exit the Test Mode immediately.
- 4) While in Test Mode, sliding the Function Switch will demonstrate the different functions of the clock:
  - a) Function Switch at OFF, it will have a continuous "Vee" sound to indicate it is in the Test Mode.
  - b) Function Switch at BUZZ, it will demonstrate the buzzer sound.
  - c) Function Switch at VIB, it will demonstrate the vibrator.
  - d) Function Switch at VIB & BUZZ, it will demonstrate the buzzer sound & the vibrator.

### ALARM CLOCK WITH BED SHAKER & USB CHARGING MODEL SBT600ss

#### Features

- |                    |                           |
|--------------------|---------------------------|
| 1) Time Set        | 8) PM indicator           |
| 2) Hour            | 9) Auto (Alarm indicator) |
| 3) Min             | 10) Tone                  |
| 4) Alarm Set / Off | 11) DC5V input            |
| 5) Alarm On/Off    | 12) Volume                |
| 6) Snooze / Dimmer | 13) USB Charging          |
| 7) Battery Cover   | 14) Vibrator Input        |

#### Setting the correct time –

- 1) Press and hold the **Time Set Button (1)** at the same time press the **Hour Button (2)** to advance the clock display rapidly until the hour is correct.
- 2) Press and hold the **Time Set Button (1)** at the same time press the **Minute Button (3)** to advance the clock display rapidly until the minute is correct.

A **lighted dot (8)** in the upper left of the clock display indicated the time is PM.

For AM time indication, the dot does not light up.

#### Setting the alarm time –

- 1) Press and hold the **Alarm Set Button (4)** at the same time press the **Hour Button (2)** to advance the alarm display rapidly until the hour is correct.
- 2) Press and hold the **Alarm Set Button (4)** at the same time press the **Minute Button (3)** to advance the alarm display rapidly until the minute is correct.
- 3) Make sure the alarm time is correctly set for AM or PM, the **lighted dot (8)** indicated PM. Please be sure to release the HOUR and MIN button before releasing the **ALARM Set button (4)**.

#### Selecting 12 or 24 hours format –

Press and hold the **Hour button (2)** for 6 seconds to enter the hour selection mode, release the **Hour button (2)** when the display shows "12 H" or "24 H". Press and release the hour button to toggle between 12 or 24 hours. If no selection is made within 4 seconds, then the last displayed hour format will be chosen and the display will return to the normal Time mode.

#### Waking to Alarm –

- 1) If only the buzzer is required, do not plug in the vibrator.
- 2) Set desired volume using the **Volume wheel (12)** for buzz sound.
- 3) If only vibrator is required, the plug in the vibrator into **Vibrator socket (14)** and turn off the volume by using the **Volume wheel (12)**.
- 4) If both the Buzzer and vibrator are required, plug in the vibrator and set the volume to the desired level.

Set alarm wake time per instructions for setting alarm time above.

#### Adjusting The Alarm Duration Time -

The default alarm duration time setting is 30 minutes, you may adjust the alarm duration between 1 minute and 59 minutes as desired. To adjust the alarm duration:

- 1) Press and hold the **AL. SET button (4)** for six seconds. The clock display will show the default alarm duration as '30', release the AL. SET button.
- 2) Press the **MIN button (3)** to adjust the alarm duration to the desired setting, from 1 minute to 59 minutes.

3) When the display is showing the desired alarm duration time, it will return to the normal clock mode after 3 seconds. The new alarm duration time setting becomes the new default and remains in the memory until you change it again.

#### Shutting Off Alarm –

Alarm can be shut off temporarily for the day or shut off permanently.

- 1) When the alarm goes off, simply press the **AL off button (4)** the alarm will shut itself off for the day and come back on the next day at the same time if the alarm setting has not changed.
- 2) To shut off the alarm permanently, slide **ON/OFF Switch (5)** to off position. The alarm will not come on until **ON/OFF Switch (5)** is set to ON position again.

#### Snooze Operation–

**SNOOZE button (6)** – press after alarm has sounded. The alarm will be turned off for 9 minutes for extra sleep. The alarm will come back again after 9 minutes. You may repeat many time up to 119 minutes .

#### Adjusting The Snooze Time

The default snooze time setting is 9 minutes but you may adjust the snooze time between 1 minute and 30 minutes as desired. To adjust the snooze time:

- 1) Press and hold the **SNOOZE button (6)** for four seconds. The clock display will show the default snooze time '9'. Release the Snooze button.
- 2) Press the **MIN button (3)** to adjust the snooze time to the desired setting, from 1 minute to 30 minutes.
- 3) When the display is showing the desired snooze time, it will return to the normal clock mode after 3 seconds. The new snooze time setting becomes the new default and remains in the memory until you change it again.

#### Volume –

To adjust the volume rotate the **Volume wheel (12)** until you have reach the desired level.

#### Tone –

Adjust the tone by rotating the **Tone wheel (10)** until you have reach the desired level.

#### Dimmer Operation

There is a 5-level brightness control for the LED display. Press the **SNOOZE button (6)** once to change the display brightness. The brightness will cycle from bright to dim then from dim to bright.

**\*NOTE\***: The dimmer function will not work, when the clock is in snooze mode.

#### Bed Shaker / Vibrator

- 1) Attach the bed shaker jack to the back of the clock at **VIBRATOR input (14)**.
- 2) Place vibrator under mattress or pillow.

#### USB Charging

This alarm clock is equipped with a convenient USB charging port on the back panel that allows you to recharge many brands of mobile phones, portable music players, and other handheld devices that are supplied with USB charging cables. The output rating of the USB charging port is 500mA.

To charge your device, connect the USB charging cable that was supplied with your handled device to the **USB charging port (13)** on the back panel of this alarm clock. Observe the charging indicator on your device to confirm that it is charging.

**Note**: Charging time will vary depending on the charging current that your device requires. Therefore some devices will take longer to charge than others.

#### Battery Backup –

To use the battery backup, you can install 1 pieces of 9V size battery at the base of the clock under the **battery cover (7)**. Please note the LED

display will not show during power failure but the clock will still be running and keeping proper time until the power is restored.


**NOTE**: For the best results it is recommended that the battery be changed every 6-8 months. If the clock is not used for a while, it is recommended to remove the battery to avoid leakage.

#### DC 5V Power Input –

Attach the DC power supply input to the back of clock labeled **DC5V --- IN (11)** .

#### GENERAL DATA

UL listed Power Adaptor: AC100-240V~, 50/60Hz  
Power Consumption: Adaptor Output: DC5V --- , 800mA  
USB Charging Port Rating: 500mA

The Lightning flash with arrowhead symbol, within an equilateral triangle, is intended to alert the user to the presence of uninsulated "dangerous voltage" within the product's enclosure that may be of sufficient magnitude to constitute a risk of electric shock to persons.

CAUTION: TO REDUCE THE RISK OF ELECTRICAL SHOCK, DO NOT REMOVE THE COVER (OR BACK) FROM THE PRODUCT. NO USER SERVICEABLE PARTS INSIDE. REFER SERVICING TO YOUR SONIC ALERT DEALER.

The exclamation point within an equilateral triangle is intended to alert the user to the presence of important operating and maintenance (servicing) instructions in the literature accompanying the appliance.

**WARNING: TO PREVENT FIRE OR SHOCK HAZARD, DO NOT EXPOSE THIS APPLIANCE TO RAIN OR MOISTURE.**

- Attention should be drawn to the environmental aspects of battery disposal.
- This icon indicates separate disposal for electrical and electronic equipment.

Normal function may be disturbed by strong electromagnetic interference. If this is the case, normal function can be resumed by simply resetting the unit following the procedure in the instructions manual. If function does not resume, try moving the unit to another location.

#### SONIC ALERT'S 2 YEAR LIMITED WARRANTY

From the moment your Sonic Alert product is purchased, Sonic Alert guarantee it for the period for 2 years. The guarantee does not cover accidents, negligence or breakage to any parts. The product is covered by the legal guarantee of conformity as provided by applicable law.

**Important**: your receipt is part of the guarantee and must be retained and produced in event of warranty claim. After sales telephone number: +44(0)1707 384438

**Geemarc Telecom, 5 Swallow Court, Swallowfields, Welwyn Garden City, Hertfordshire, AL7 1SB, UK**

#### RECYCLING DIRECTIVES

The WEEE (Waste Electrical and Electronic Equipment) has been put in place for the products in end of life are recycled in the best way.

When this product is out of order, please do not put it in your domestic waste bin. Please use one of the following disposal options:

- Remove the batteries and deposit them in an appropriate WEEE skip. Deposit the product in an appropriate WEEE skip.
  - Or hand the old product to the retailer. If you purchase a new one, they should accept it.
- Thus if you respect these instructions you ensure human health and environmental protection.

**CE Declaration**: Hereby Sonic Alert declares that this product is in compliance with the essential requirements and other relevant provisions of the Electromagnetic compatibility (EMC) 2004/108/EC, and of Low Voltage Directive (LVD) 2006/95/EC (if applicable).

**The declaration of conformity may be consulted at [www.geemarc.com](http://www.geemarc.com)**

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