

TESTIMONIALS



“The mCareWatch saved my mother’s life.”

Lawrence Hawthorne | QLD

“It’s not very often that a product changes people’s lives, but the mCareWatch certainly does.”

Paul Neville | NSW

“Mum had a fall at home and was able to press the SOS button on the watch which automatically called me and I was able to get her help.”

Tim Graham | NSW



NEW

The first 4G SIM-enabled, water resistant smart watch, personal emergency alarm, GPS tracker and mobile phone providing unrivalled support for **caring, safety and wellbeing**



mCareWatch HW2

AFFORDABLE SAFETY, INDEPENDENCE + PEACE OF MIND

Achieve new levels of caring, safety, peace of mind and functionality with the mCareWatch HW2, which we can preconfigure to connect wearers, carers and their devices to the mCareWatch ConnectiveCARE platform.

The mCareWatch HW2 is an all-in-one watch, mobile phone, personal emergency alarm, GPS & activity tracker. It is a medical alert, can receive customised medication and appointment reminders, can be linked to 24-hour emergency personal monitoring services or to auto-dial a family or friend in case of emergency.

With 4G network connectivity, the mCareWatch HW2 is the most advanced, easiest to use, program and recharge wearable device on the market.

It's the ideal assistive technology to support the safety, health and wellbeing of people who have:

- Dementia/Alzheimer's disease
- Disabilities such as autism, epilepsy or developmental disorders
- Heart disease, diabetes or other chronic conditions that cause concern about wellbeing
- Recent hospital discharge or transitioning from hospital to home



FEATURES

TIME



Tells the time, day and date with analogue & digital watch faces available

SOS CALLS



SOS button that dials pre-set contact number or 24/7 monitoring

LOCATION FINDER



Locate the watch wearer via GPS, WiFi and GSM (mobile towers)

CALLS



Watch has it's own mobile number which can make and receive calls

WATER-RESISTANT



Watch is water resistant and can be worn in the shower



Setup a safe perimeter zone and be notified if the wearer exits this safe zone

MOBILE APP & WEB PLATFORM



Carers use the mobile app and online platform to locate watch, update SOS contact, setup reminders, etc.

NON-MOVEMENT



Carers receive an alert when wearer has not had any movements during specific times during the day

CUSTOMISED REMINDERS



Carers can setup customised medication and appointment reminders for the wearer

FAMILY SPEED DIAL



Wearer can call up to 7 contacts by simply pressing photo of contact

CAPTURE BLOOD PRESSURE



Watch connects via Bluetooth with separate biometric devices to measure blood pressure and weight

WELLBEING CHECK-IN



Enables the watch to check in on the wearer by asking a series of questions to ensure they are safe