

Timer Comparisons

| | Approx. Vol.* | Product | Alarms | Compartments | Alarm Duration | Model | Price |
|---|---|---|---|--------------|--|--------------------------------------|--------------------|
|  | 35-40 db + (Conversation, Radio) +Vibrating Alarms (suitable for people with hearing difficulties) | VibraLITE Watches | Up to 2, 8 or 12 p/day depends on model | 0 | 20 seconds unless stopped sooner | Various | \$85.95 - \$229.95 |
|  | 50-70 db (TV, Vacuum Cleaner, Car 45klm-25ft) A: 50-60 db B: 55-60 db A+B: 65-70 db | Automated Pill Dispenser (Solid White Lid or Clear Lid) | Up to 4 p/day | 28 | Alarm and blinking light for up to 30 minutes. Pills stay in window until next alarm due | TabTimer TT4-28SW or TT4-28CL | \$149.95 |
|  | 65-70db (TV, Vacuum Cleaner, Car 45klm-25ft) | Pill Reminder 37 | Up to 37 p/day | 7 | 10 seconds then every 10 seconds for 4 minutes. Then 'Missed Pill' message will flash on LCD display | TabTimer TT37-7 | \$49.95 |
|  | Hi/Lo 45-55db (Average Home, Office) | TabTime Super 8 | Up to 8 p/day | 8 | 40 seconds then LED flashes until opened | TabTimer TTS8-1 | \$49.95 |
|  | 45-50 db (Average Home, Office) | Pill Box Reminder | Up to 4 p/day | 3 | 3 times then 3 times every minute until confirmed | TabTimer TT4-3 | \$26.95 |
|  | 40-45 db (Average Home) | Reminding Clock | Up to 4 p/day | 0 | 3 times then 3 times every minute until confirmed | TabTimer TT4-0SQ | \$31.95 |
|  | 30-40 db (Quiet Street Traffic, Bird Calls) NOT suitable for people with hearing difficulties | Anywhere mini p-Clock | Up to 4 p/day | 0 | 3 times then 3 times every minute until confirmed | TabTimer TT4-0RND | \$19.95 |
|  | 30-40 db (Quiet Street Traffic, Bird Calls) NOT suitable for people with hearing difficulties | Pill Box Alarm 4 | Up to 4 p/day | 2 | 60 seconds until cancelled | TabTimer TT4-2Y / 2C | \$22.95 |

* Approximate volumes established at 1m unobstructed distance using digitech micro sound level meter comparison examples provided as a guide only.