

PRODUCT INFORMATION



Stay on-time with TabTimer Reminders!

Since 2009 remembering to administer medicines on-time and carry out personal care tasks at the right time has been made easier using assistive technology and assisted listening reminder devices available from exciting Australian reminders company - TabTimer™.

TabTimer™ supplies a handy range of quality, innovative and affordable reminder devices which can be easily set to alarm when medicines, vitamins, or other personal tasks are due – serving as a handy reminder to help patients, their relatives, carers or nursing staff to achieve medication compliance or remember that easy forgotten personal care task.

TabTimer's range includes pill box timers with multiple alarms, electronic pill boxes, lockable automated pill dispensers, handy pocket reminders, fridge magnet reminders, vibrating or talking alarm watches and clocks. Large displays and controls make it easy for the elderly, disabled or vision impaired person to read and operate. A selection of TabTimer's devices can also be set to flash or vibrate only - making the process of remembering to take medicines or perform daily living tasks discrete and timely.

For more information visit www.TabTimer.com.au or call 1300 TAB TIMER (1300 822 846) in NZ (0800 822 846)

TabTimer™ helps keep medicines on time™



- ENDS -

Publicity Contact:-

Jonathan Allott
Managing Director, TabTimer Pty Ltd
P - AU: 1300 TAB TIMER (1300 822 846)
p - NZ: 0800 TAB TIMER (0800 822 846)
F - +61 1300 883 749
E - jonathan@tabtimer.com.au

Information Sheet

DVA - Rehabilitation Appliances Program (RAP)

If you are a veteran who is eligible for support from the Department of Veterans' Affairs and hold a gold or white card, you may be able to claim one of our products under the Rehabilitation Appliances Program (RAP). Our DVA Supplier Number is ; S0028918 .

To access DVA funding, you will need to speak to your local medical officer, GP, or health professional (e.g. an occupational therapist, specialist) to establish if you are eligible and have them complete the necessary applications.

TabTimer Products on DVA Contract

6 of our products are now on contract with DVA and available via Allianz, Bright Sky and Country Care. Occupational Therapists and other allied health staff can order the contracted products from the suppliers' websites.

The products currently available on contract are the; Automated Pill Dispenser, Vibe-5, VibraLITE MINI, MedCenter Talking Clock, and MedCenter Your Minder Recordable Clock.

If you need another product, this is possible but it is obtained via a Prior Approval from DVA, using a Direct Order Form. Additional details on that process are also shown below.

Prior Approval Process

What is RAP?

DVA – (RAP) assists eligible persons to be as independent and self-reliant as possible in their own homes. The provision of aids and appliances is intended to minimise the impact of disabilities and maximise quality of life and independence. RAP provides aids and appliances according to assessed clinical need. Eligible persons are assessed by an appropriately qualified health professional, usually an occupational therapist, GP, physiotherapist or registered nurse. Where appropriate, the health professional will also provide instruction in the use of the aid/appliance. Items issued under RAP are generally designed specifically for people with an illness or disability.

Accessing RAP

Where an aid or appliance is assessed as clinically necessary, the assessing health professional must complete a Direct Order Form (D0992). Forward the form, along with a clinical justification, to TabTimer so we can prepare a quote and we will then forward the application to DVA Rehabilitation Appliance Program (RAP), who will then consider the application and, if approved, will forward the order directly to TabTimer.

Non-Schedule items







TabTimer products that are not listed as an available item on the DVA Schedule of Equipment can be obtained via a special application for a 'non-Schedule item'. One-off requests for non-Schedule items may be considered where an assessed clinical need cannot be adequately addressed by any current Schedule item.

If you require a quote, you can contact TabTimer and request one, or send us your completed Direct Order Form -D0992 along with a clinical justification, and we will submit it to DVA for you with a copy of our quotation.

For more information contact TabTimer.

DVA Contracted Items

The following TabTimer products are now available via funding from the Department of Veterans' Affairs. The items can be prescribed by; Local Medical Officers or GPs, Occupational Therapists, Physiotherapists, Registered Nurses, and Medical Specialists via an order for Veterans from the suppliers websites shown in the table below.

	<p>DVA FUNDING</p> <p>Automated Pill Dispenser – TT6-28SC Up to 6 alarms per day - 28 storage compartments <i>Load up your medicine for 7 to 28 days</i> <i>Solid White Lid or Clear Interchangeable Lids included in pack</i></p> <ul style="list-style-type: none"> • See further details on page 9 <p>Order now from; Allianz Bright Sky Country Care</p> <p style="text-align: right;">RRP AU \$199</p>
	<p>DVA FUNDING</p> <p>VibraLITE 12 Mini - Vibrating Reminder Watch Up to 12 Daily Alarm Settings <i>Smaller Case & Watch Band - For Small Wrists of All Ages</i></p> <ul style="list-style-type: none"> • See further details on page 23 <p>Order now from; Allianz</p> <p style="text-align: right;">RRP AU \$109</p>
 <p>TTC-MCTALK AU \$ 85</p>  <p>Optional adapter MP3310 AU \$ 29</p>	<p>DVA FUNDING</p> <p>MedCenter Talking Alarm Clock -TTC-MCTALK Talking Alarm Clock and Medication Reminder <i>Up to 4 medication reminders per day</i></p> <ul style="list-style-type: none"> • See further details on page 12 <p>Order now from; Allianz Bright Sky</p> <p style="text-align: right;">MCTALK RRP AU \$85</p> <p style="text-align: right;">ADAPTOR RRP AU \$29 Allianz</p>
 <p>TTC-MCREC AU \$99</p>  <p>Optional adapter MP3310 AU \$ 29</p>	<p>DVA FUNDING</p> <p>Your.Minder Personal Alarm Clock –TTC-MCREC Talking Alarm Clock with Six Personal Messages <i>Up to 6 reminders per day</i></p> <ul style="list-style-type: none"> • See further details on page 12 <p>Order now from; Allianz Bright Sky</p> <p style="text-align: right;">MCRRC RRP AU \$99</p> <p style="text-align: right;">ADAPTOR RRP AU \$29 Allianz</p>

DVA Contracted Items



Amplicall 1 with 50 Bundle - TT-AC1-50B
Wireless door bell with Vibrating Ringing Indicator Alerter
for deaf and hearing impaired

TT-AC1-50B – is replaced by the door bell
TT-VB2+1 , awaiting DVA/Allianz to update this item

- See further details on page 33

Order now from;



RRP AU \$169

Allianz



MemRabel 2 (v3.1) - TTC-MEMRABEL2
Audio Visual Orientation/Calendar Alarm Clock
Perfect for people with Alzheimer's, dementia, memory
impairment, cognitive impairment, low vision, or everyday
forgetfulness

- See further details on page 18

Order now from;

DVA
FUNDING

Allianz

RRP AU \$299



VISO10 - TTC-VISO10
Large Easy Read Alzheimer's Dementia Memory Clock
Non Abbreviated Display of Day, Time, Date, Month and Year
Perfect for people with memory impairment

- See further details on page 17

Order now from;

DVA
FUNDING

Allianz

RRP AU \$189

Other TabTimer items

**For any other TabTimer product please
complete a Direct Order Form to obtain prior
approval as detailed on page 2 .**

DVA
FUNDING

Timer Comparisons

	Approx. Vol. *	Product	Alarms	Compartments	Alarm Duration	Model	Price
	35-40 db +vibes (Conversation, Radio) +Vibrating Alarms (suitable for deaf people & with hearing difficulties)	Vibrating Watches	Up to 8 or 12 p/day depends on model	0	20 seconds unless stopped sooner	Various	AU \$109 - \$199
	35-40 db (Conversation, Radio) Sound and/or Vibration alarms	CADEX VibraPLUS 8 Watches	Up to 8 p/day	0	20 seconds unless stopped sooner	Various	AU \$179 - \$199
	50-70 db (TV, Vacuum Cleaner, Car 45km-25ft) A: 50-60 db B: 55-60 db A+B: 65-70 db	Automated Pill Dispenser (Solid White Lid and Clear Lid)	Up to 6 p/day	28	Alarm and blinking light for up to 30 minutes. Pills stay in window until next alarm due	TabTimer TT6-28SC	AU \$199
	Hi / Lo 55-75db (Average Home, Office)	TabTime Super 8	Up to 8 p/day	8	40 seconds then LED flashes until opened	TabTimer TTS8-8	AU \$59
	40-55 db +vibes (Average Home) +Vibrating Alarms (suitable for deaf people & with hearing difficulties)	Vibrating 12 Alarm pager	Up to 12 p/day	0	20 seconds plus one 5 min snooze per alarm possible by pressing Lite	TabTimer TT12-0SQ	AU \$119
	45-50 db (Average Home, Office)	Pill Box Reminder	Up to 4 p/day	3	3 times then 3 times every minute until confirmed	TabTimer TT4-3	AU \$34
	40-45 db (Average Home)	Reminding Clock	Up to 4 p/day	0	3 times then 3 times every minute until confirmed	TabTimer TT4-0SQ	AU \$34

* Approximate volumes established at 1m unobstructed distance using digitech micro sound level meter. Comparison examples provided as a guide only.

TabTimer Product Listing

Timers



Reminding Clock - TabTimer TT4-0SQ

Up to 4 alarms per day - square reminder

Aural or visual alert with magnetised clip

- Place on the fridge or other convenient location for medicine reminder
- Up-to 4 alarms per day - Easy to set and use
- Alarm sounds 3 times then every minute until confirmed
- 12 hour clock with 'PM' indicator for PM LCD display
- 2 x "AAA" size batteries included
- Unit Dimensions: 69 X 49 X 15mm , vol. approx. 40-45db
- 12 Month Limited Warranty by product source

RRP AU \$34



TabTimer 8 Alarm Timer - TabTimer TT8-0SQ

Up to 8 alarms per day - square reminder

Aural or visual alert with magnetised clip

- Place on the fridge or other convenient location for medicine or other reminder
- Up-to 8 alarms per day - Easy to set and use
- No need to reprogramme, the set alarms will repeat daily
- Countdown timer for regular doses (up to 99 hours and 59 mins)
- Alarm sounds for 60 seconds when reminder time arrives
- 3 large buttons control all functions, Large Easy read LCD
- 12 hour clock with 'PM' indicator for PM or
- Clock with 24 hour display function (military style)
- Magnet on back, stand on desk, or hang on wall
- 1 x "AAA" size battery included
- Unit Dimensions: 80 X 56 X 20mm , vol. approx. 40-45db
- 12 Month Limited Warranty by product source

RRP AU \$49



Vibrating 12 Alarm pager - TT12-0SQ

Up to 12 alarms per day - vibrating (or sound) reminder

Fits in your Pocket, Purse, or use included Belt Clip

- STRONG Vibration Alarm (or Sound Alarm) - Up to 12 Daily
- Fits in your Pocket, Purse or use included Belt Clip
- Count down timer and more! (up to 999 hours)
- Some initial programming required
- No Monthly FEES or Paging Service required
- Alarm sounds for 20 seconds with one 5 min snooze per alarm
- Unit Dimensions: 55 X 40 X 15mm, vol. approx. 40-45db
- 12 Month Limited Warranty by source (excludes battery)

RRP AU \$119

Pill Box Alarms

Pill Box Reminder - TabTimer TT4-3

Up to 4 alarms per day - 3 storage compartments

Pocket sized case, good for carrying



- Up-to 4 alarms per day
- 3 compartments for storing different medicines
- Easy to set and use
- When due, alarm sounds 3 times then every minute until stop pressed
- 12 hour clock with 'P' indicator for PM
- LCD display, 1 x "LR44" size battery included
- FDA grade polypropylene
- Compartment Dimensions : 40 x 25 x 20mm
- Total Unit Dimensions: 120 X 45 X 23mm, vol. approx. 45-50db
- 12 Month Limited Warranty by product source
- Made in Taiwan

RRP AU \$34

TabTime Super 8 - TabTimer TTS8-8TB

Up to 8 alarms per day - 8 storage compartments

One of the most effective and efficient daily pill timer and medicine reminder on the world market.



**DVA
FUNDING**

- Can be set at exact times when reminders are due
- Countdown timer for regular doses (up to 99 hours and 59 mins)
- Alarm sounds for 40 seconds when reminder time arrives
- LED continues to flash until opened indicating 'missed dose'
- Trigger switch stops alarm/flashing once device is opened Easy open magnetic catch for weak or shaky hands
- Adjustable 'Hi/Lo' Volume (Volume level approx 45-55 db)
- 3 extra large control buttons, Large display screen
- Easy access battery compartment for easy replacement
- Clock with 24 hour display function (military style)
- 1 x "CR2032" Size battery included, battery indicator
- Compartment Dimensions : 25 x 25 x 13mm
- Total Unit Dimensions: 110 x 60 x 30 mm, Made in China

RRP AU \$59

Now, it's easy to simplify complex dosing schedules using this proven product - the VitaCarry Pill Reminder. Always know when and what you need to take because it automatically reminds you. Audible and visual alarms help you keep track throughout your busy day. Programmable in seconds, this handy device can remember up to 4 or 7 daily alarms and as an assistant - even let you know if you've missed a dose. Slips easily into a handbag or glove box, you will also appreciate the sturdy, food grade plastics that keep expensive contents safe. Secure, easy-to-use latch provides easy access to the contoured compartments. Requires two AAA batteries (included). View critical information at a glance through clear, scratch resistant lens. Always know when your next dose is scheduled (alarm time dominates the screen). Improves medication adherence and quality of life. Program up to four or seven daily alarms. Audio alarm can be turned off for more discreet visual alarm feature.



TT-VC4-4 RRP AU \$59
VitaCarry 4 Compartment Pill Box
with 4 Alarm Timer



TT-VC4-0 RRP AU \$49
VitaCarry 4 Alarm Timer Only



TT-VC0-5 RRP AU \$19
VitaCarry 5 Compartments Only

VitaCarry Gasketed 4 Compartment Pill Box (White) with 4 Alarm Timer
4 Dose Large Capacity Organizer and Medication Reminder
Load up your medicine for 4 doses with 4 daily alarms

- 4 Compartment Pill Reminder with 4 times per day alarm timer
- Audible beeping alarm, vibration, and flashing display
- Gasketed to keep air tight and protect your medicines against damage from air and moisture
- Alarm Timer can easily be removed to make it dishwasher safe
- Made of sturdy food grade plastic for maximum protection
- Contoured compartments for easy pill removal
- Missed dose indicator, Easy to use
- Dimensions 15cm x 10cm x 3cm
- Batteries: 2 x AAA included
- 12 Month Limited Warranty replacement by product source

4 Compartment Pill Box with 4 Alarm Timer RRP AU \$59
4 Alarm Timer Only - RRP AU \$49
5 Compartment Pill Box Container Only - RRP AU \$19

Automated Pill Dispenser

Automated Pill Dispenser – TT6-28SC

Up to 6 alarms per day - 28 storage compartments

Load up your medicine for 7 to 28 days

Solid White Lid or Clear Interchangeable Lids included in pack



- Automated pill dispenser
- 28 Pill compartment design -Up-to 6 alarms per day
- Alarms for up to 1 – 2 – 3 – 4 – 5 or 6 times per day
- Compartments can be loaded with medicine for 7 to 28 days, depending on number of doses per day
- Lockable pill dispenser, ensures only the right dose is accessible
- Simple to set up and easy to use, with long duration alarms
- Loud alarm and blinking light for up to 30 minutes
- 3 alarm tones plus silence to accommodate different hearing levels
- Alarm automatically resets when medicine taken
- Medicines not taken remain in window until next alarm
- Clock with 12 or 24 Hour Display Function
- Alarms can be set at any time during the day – up-to 6 per day
- Locking key system, Built in test function, Low battery indicator
- 4 x "AAA" Size batteries included , vol. approx. 50-70db
- Pharmacies can use tamper evident labels if required, to meet QCPP - Quality Care Pharmacy Program requirements for a dose administration aid. 4x Tamper Evident labels supplied.
- Unit Dimensions: Diameter 192.5mm x 52mm High
- 12 Month Limited Warranty replacement by product source

RRP AU \$199

**DVA
FUNDING**

Order now from;

[Allianz](#)

[Bright Sky](#)

[Country Care](#)



Spare Tray Set for Automated Dispenser – TT6-28XT

- Additional trays are also available
Some customers leave a tray at the pharmacy while using the other to save time when medicines need refilling

RRP AU \$35



**AUTOMATED
PILL DISPENSER**
with interchangeable lids
(clear lid & solid lid)

**PLUS
SPARE TRAY
BUNDLE**

with spare tray, dust cover
lid, dosing ring templates,
and spare key

Automated Dispenser and Spare Tray Set Bundle - TT6-28SCB

- Contains Automated Dispenser and spare tray set

RRP AU \$219

Medready Automatic Pill Dispenser

MedReady MR-1700FL Automatic Pill Dispenser

Up to 4 alarms per day - 28 storage compartments

Load up your medicine for 7 to 28 days



- Fully automatic controlled pill dispenser – WITH SOLID LID
- 28 Pill compartment design (29 compartments with 28 usable)
- Up-to 4 alarms per day, 1 – 2 – 3 or 4 times per day
- Compartments can be loaded with medication for 7 to 28 days, depending on number of doses per day
- Lockable pill dispenser, tamper resistant, ensures only the right dose is accessible at the right time
- Large 25mm LED displays time and aids setting alarms
- Simple to set up - as easy as a digital clock
- Loud & long duration alarms up to 250 minutes (3 - 250 minutes)
- When alarm time is reached device beeps 5 seconds on, 1 second off until medicine door is opened
- Optional beep/chirp every minute if dose is missed until dose taken or next dose alarm is reached
- Powerful alarm, 90dB, 400Hz and Flashing Light.
- Audible and visual alarms can be independently turned on/off
- Alarm automatically resets when medication taken (no need to program daily unless change required)
- Medicines not taken will remain in window until next alarm
- Missed doses will move around and become inaccessible upon next alarm, allowing monitoring of adherence to regimen
- Clock with AM/PM Display Function, Built in test function
- AC adaptor and Rechargeable Battery included
- Unit Dimensions: Diameter 270mm x 55mm High
- Each compartment is approx. 35mm x 20mm x 10mm (large enough for approx. 5-6 large vitamin tablets or approx. 10-12 panadol sized tablets.)
- Pill compartments are sloped for easy removal of pills, divisions keep medicines in place even if device tipped
- Card holder on the bottom of the unit for a printed copy of the patient's prescriptions for all medicines in the dispenser.
- 12 Month Limited Warranty replacement by product source (excludes battery), Made in USA

RRP AU \$ 299

MedReady 1700PRN Automatic PRN Pill Dispenser

PRN dispensing simply limits the time between doses

28 storage compartments - Strong locking lid for added security



- Dispenses medication for the typical PRN medication regimen
- PRN dispensing simply limits the time between doses
- This dispenser lets you preset the minimum time between doses from 15 minutes to 24 hours
- Upon accessing a single dose of medication, a visible preset countdown timer starts which indicates how long before the next dose is available
- When this timer reaches 0:00, the next available dose becomes available
- Note that to make it a true PRN (as needed) device, we do not sound any alarms indicating the next dose is available
- 12 Month Limited Warranty replacement by product source (excludes battery), Made in USA

RRP AU \$ 299

Medready – Internet connected for monitoring & alerts



MedReady MR-357FL Automatic Pill Dispenser with built in 3G SIM card and Flashing Light (Capable of SMS & Email Alerts*)

Up to 4 alarms per day - 28 storage compartments

Load up your medicine for 7 to 28 days

Caregivers can change settings over the Internet and receive alerts if pills missed*

*Includes built in 3G SIM for Remote Self Monitoring, SMS & Email Alert Notifications when scheduled dose has not been taken. *Requires monthly subscription fee (\$19.95 p/mth see MedReady Subscription MR-SUB). No long term contract is required and you may cancel at any time.



Flashing light for the hearing impaired



- Fully automatic controlled pill dispenser - WITH SOLID LID
- 28 Pill compartment design (29 compartments with 28 usable)
- Up-to 4 alarms per day, 1 – 2 – 3 or 4 times per day
- Compartments can be loaded with medication for 7 to 28 days, depending on number of doses per day
- Lockable pill dispenser, tamper resistant, ensures only the right dose is accessible at the right time
- Large 25mm LED displays time and aids setting alarms
- Simple to set up - as easy as a digital clock
- Loud & long duration alarms up to 250 minutes (3 - 250 minutes)
- When alarm time is reached device beeps 5 seconds on, 1 second off until medicine door is opened
- Optional beep/chirp every minute if dose is missed until dose taken or next dose alarm is reached
- Powerful alarm, 90dB, 400Hz and Flashing Light.
- Audible and visual alarms can be independently turned on/off
- Alarm automatically resets when medication taken (no need to program daily unless change required)
- Medicines not taken will remain in window until next alarm
- Missed doses will move around and become inaccessible upon next alarm, allowing monitoring of adherence to regimen
- Clock with AM/PM Display Function, Built in test function
- AC adaptor and Rechargeable Battery included
- Unit Dimensions: Diameter 270mm x 55mm High
- Each compartment is approx. 35mm x 20mm x 10mm (large enough for approx. 5-6 large vitamin tablets or approx. 10-12 panadol sized tablets.)
- Pill compartments are sloped for easy removal of pills, divisions keep medicines in place even if device tipped
- Card holder on the bottom of the unit for a printed copy of the patient's prescriptions for all medicines in the dispenser.
- 12 Month Limited Warranty replacement by product source (excludes battery), Made in USA

RRP AU \$499

Remote Self Monitoring

- Caregiver access to control device remotely via the Internet
- Ability to program and self monitor device remotely
- Create Alert Notifications via SMS and/or email if scheduled dose is missed, power is disconnected or an error occurs
- View status of device, log and remaining doses left in machine
- Change alarm times, alarm duration, turn alarms and flashing light on/off, on website
- Remotely test the device and provide an early dose if required
- View detailed reports

RRP AU \$18 - \$22 p/mth

Talking Clocks

"Good morning!
The time is
10.30 AM,
Friday, January

"Good morning!
Please take your
morning pills
for the 20th.
BEEP BEEP BEEP"



The MedCenter System – TTC-MCSYS Monthly Organizer and Talking Pill Reminder

Load up your medicine for 31 days

- Set up to 4 ALARMS DAILY to alert you when your medication is due
- TALKING ALARM alerts users when it's time to take medication
- A repeating reminder notifies the time, date and dose
- "Good morning please take your morning pills for the 20th"
- Press Alarm Acknowledged to confirm
- Talking Alarm Clock tells you when next dose is due
- TRI-DATE VERIFICATION: matching the visual dates on pill boxes and clock with audible alerts assures accuracy
- LOUD & EXTRA LOUD sound settings ensure the alarm is heard
- RED & GREEN COLOR CODING shows when daily dose is complete
- Dimensions;
Base- 250mm W x 180mm H x 200mm D
Clock – 125mm W x 65mm H x 45mm D
Individual Pill compartment – 25mm x 25mm x 20mm H
Daily Pill Box (4 compartments) – 100mm x 25mm x 25mm

RRP AU \$159



TTC-MCTALK AU \$ 85



Optional adapter MP3310 AU \$ 29

[Allianz](#)

MedCenter Talking Alarm Clock -TTC-MCTALK Talking Alarm Clock and Medication Reminder

Up to 4 medication reminders per day

- TALKING ALARM alerts users when it's medicine time
- Set up to 4 ALARMS DAILY to alert you when your medication is due
- LOUD & EXTRA LOUD sound settings ensure alarm is heard
- EASY SET controls make setting the clock even easier
- Dimensions; 125mm W x 65mm H x 45mm D

RRP AU \$85

**DVA
FUNDING**

Order now from;

[Allianz](#)

[Bright Sky](#)



TTC-MCREC AU \$99



Optional adapter MP3310 AU \$ 29

[Allianz](#)

Your.Minder Personal Alarm Clock –TTC-MCREC Talking Alarm Clock with Six Personal Messages

Up to 6 reminders per day

- RECORD UP TO (6) OF YOUR OWN ALARM MESSAGES to alert users when it's time to take medication, get out of bed, eat, go to work, walk the dog, etc.
- Set up to (6) DIFFERENT DAILY ALARM times
- LOUD speaker volume ensures the alarm is heard
- Powered by (2) AA Batteries (included) or option for AC Adapter purchased separately
- Dimensions; 125mm W x 65mm H x 45mm D

RRP AU \$99

**DVA
FUNDING**


Order now from;

[Allianz](#)

[Bright Sky](#)

Talking Clocks

Reminder Rosie))) PERSONAL VOICE REMINDER



"I love you, Mom. It's time to take 2 red & 1 white pill"
8:30 AM every day



"Dad, your doctor appointment is in one hour"
1:00 pm, Fri Dec. 4, 2015

"I just wanted to let you know that after six months addition to my parents daily lives of use, Reminder Rosie has been very successful. It has enabled them to maintain the schedule of the assisted living residence and, I think, my parents seem to enjoy hearing me and my sister's voice giving them the reminders. Thanks for bringing such a helpful device to market at a reasonable price Rosie has been a great help."

— CN, Smithville, Ontario

Features

Hands-Free Trigger

Reminder Rosie is exceptionally easy to use! Many devices require fiddling with buttons and switches to operate. With Reminder Rosie, you just start talking and Reminder Rosie is ready!

Simple Yet Sophisticated

Behind a simple & familiar clock interface is a robust technology. Reminder Rosie processes spoken language almost as easily as humans do. By carefully listening for words or phrases, Rosie is ready 24/7 to execute your requests.

Custom Reminders

Reminders can be programmed in any language to provide greater comfort, familiarity, and to improve task completion.

Multiple Reminder - Any Time or Date

With Rosie, you can record up to 25 reminders for everyday, weekly or ANY date. Review the day's reminders with a simple voice command so you stay in control.

Emergency Back-up System

Rosie will work even during power outages. Rosie will maintain its daily reminders and all other critical features.

Reminder Rosie – TTC-ROSIE-AU

Personal Voice Activated Reminder Alarm Clock

Record up to 25 reminders for everyday, weekly or any day

Reminder Rosie is an advanced, voice controlled reminder system that solves the very real daily challenges of memory loss. Inspired by the creator's father who suffered from dementia, Rosie is a personalized, voice controlled reminder system which is so simple that older adults love to use it! With its senior-friendly clock interface, Reminder Rosie is a simple, hands-free, inexpensive solution to remember medication, appointments, and every-day tasks. Never forget a task again!

Benefits

Maximize Independence & Engagement

- Prolongs senior independent living and minimizes dramatic life changes
- Helps seniors maintain their healthy routines
- Assists healthcare providers and family members to keep track of older adults schedule and required tasks
- Requires no behavioural change to ensure easy adoption and use by Seniors

Relieve Caregiver Stress

- Minimizes or even eliminates the "nag factor"
- Reminder Rosie eases stress by taking over the mundane tasks like reminders from daily responsibilities so caregivers can focus on what matters most – extending the best quality of care to the older adult

Assist Visually Impaired

- Using state-of-the-art voice activated technology, Reminder Rosie can operate almost entirely by voice, thus greatly reducing the need for the older adult to fiddle with buttons

Increase and Sustain Medication Adherence

- Reminds adults to take their medications in a customized manner throughout the day
- Helps improve adherence over a sustained period of time
- Minimizes medication disruptions
- Reduces hospitalizations that are linked to medication issues

Provide Personalized Care

- Creates needed custom care for the older adult with minimal oversight
- The platform is voice customized for greater adherence

Specifications

- Hands Free Voice Activated Reminder Aid
- Revolutionary Speech Recognition
- Voice or manual operation
- Record 25 reminders, 6 seconds in any language
- Specific Time and Date, Daily, Weekly, Monthly & Annually
- Daily Calendar Organizer Digital Alarm Clock
- A Simple Solution to Organize your day- Remember Meds, Appointments, Household & Personal Tasks, Birthdays & Much More.
- Large 2 inch LED Digits visible from any angle
- Hi/low volume and dim settings; Alarm on/off switch
- Doubles as a voice alarm clock
- User says "reminder off" or press down to turn off
- Understands 25 spoken commands in conversational English
- Snooze function 9 minutes
- Works offline
- Power: AC 240 volts, 60Hz, AU plug
- Battery back-up: 3AAA not included
- Clock Dimensions/Weight: 20W x 10.5H x 5D
- Warranty: 1-year limited

RRP AU \$169

Memory Aid Clocks



8" Digital Calendar Day Clock - TTC-DC8001 Orientation Clock for Alzheimer's & Memory Impairment Non Abbreviated Display of Day, Time, Date, Month and Year Perfect for people with memory impairment, New Model - Has up to 5 basic Alarms per day and up to 3 Medicine Reminders

Reduce the workload of carers

The Dementia Day Clock is aimed at helping people living with Dementia by simply and clearly displaying the time of day as morning, afternoon, evening or night. People living with Alzheimer's or Dementia can lose the ability to recognise what the time is and even whether it is day or night. This can cause many challenges for them and for their primary carers. By using this digital calendar orientation day clock it can reduce the challenges.

Help people with dementia live in a daily routine

The aim of the day clock is to help people with dementia to be able to link the time of day to their daily routine. The Day Clocks' clear and simple display has proven to be helpful for people with dementia or Alzheimer's

Application Scenario

The Day Clock can be used in the living room, the bedroom, the study, the kitchen and classroom. It can be standing or wall mounted. The high gloss framing is good looking to suit your home environment perfectly.

9 Languages

The time display and menus can be displayed in one of 8 available languages; English, French, German, Dutch, Italian, Polish, Spanish Portuguese and Finnish

Time Mode / Date Mode

The time can be displayed in 12 hours or 24 hours format. The date can be displayed as Month-Date-Year or Date-Month-Year

User Interface Feature

The clock can be used in the standard mode displaying the Day, Time, Date, Month and Year. Alternatively, you can choose the analogue clock setting, or a coloured background. Press and hold the play button on the back to switch to the desired display.

Alarms

The latest model has; • up to 5 basic alarms (Monday to Friday, Saturday to Sunday, or Monday to Sunday) and • up to 3 Medicine Reminders (one time only, or repeating daily) The user can stop the alarm by pressing any button

Automatic Brightness

The clock has the function to automatically adjust the display depending on the time of day. High, Medium, Low or Auto. In Auto, the display brightens during Day Mode after 7am (50cd/m²) and automatically dims during Night Mode after 7pm (10cd/m²).

Features:

- 8 inch LED Backlight Digital Display Clock
- Day, Time, Date, Month and Year displayed in full wording - non abbreviated, Optional display in analogue clock or coloured background
- Automatic brightness adjustment setting; Day Mode after 7am (50cd/m²) and Night Mode after 7pm (10cd/m²)
- Adjustable volume (levels 00-30)
- New Model - has up to 5 basic alarms (Monday to Friday, Saturday to Sunday, or Monday to Sunday)
- Up to 3 Medicine Reminders (one time, or repeating daily)
- The user can press any button to stop the alarm
- Display in one of 9 Languages; English, French, German, Dutch, Italian, Polish, Spanish, Portuguese, or Finnish
- 12 hour or 24 hour time format
- Date display format; Date-Month-Year or Month-Date-Year
- Time Display; Morning 6am-11am, Afternoon 12pm-5pm, Evening 6pm-11pm, Night 12am-5am
- Resolution 800 x 600 high definition display
- Display Ratio 4:3
- No glare within 15 feet
- Dimensions 217mm W x 172mm H x 25mm D
- Approximate dimensions of the characters; °Day characters - 25mm H x 12.5-18mm W
- Time numerals - 38mm H x 20mm W
- Date characters - 10mm H x 8mm W
- Weight 410g
- Accuracy +/- 30-60 secs p/mth
- Electricity operated only. Internal battery retains the clock time in the event of a power failure. However, the clock must be plugged in and powered by electricity to operate.
- 5V/1A Power Adaptor supplied



RRP AU \$149



8" Digital Wi-fi Weather Station Orientation Clock - TTC-DC8001-W wifi

Orientation Clock for Alzheimer's & Memory Impairment
Non Abbreviated Display of Day, Time, Date, Month and Year
Perfect for people with memory impairment, Up to 5 basic Alarms p/day & up to 3 Medicine Reminders , plus local weather information

Reduce the workload of carers

The Dementia Day Clock is aimed at helping people living with Dementia by simply and clearly displaying the time of day as morning, afternoon, evening or night. People living with Alzheimer's or Dementia can lose the ability to recognise what the time is and even whether it is day or night. This can cause many challenges for them and for their primary carers. By using this digital calendar orientation day clock it can reduce the challenges and let the person know what the weather is like outside so they can dress appropriately.

Help people with dementia live in a daily routine

The aim of the day clock is to help people with dementia to be able to link the time of day to their daily routine. The Day Clocks' clear and simple display has proven to be helpful for people with dementia or Alzheimer's. When connected to wi-fi they can have up to date and forecast weather information , temperature and humidity.

Application Scenario

The Day Clock can be used in the living room, the bedroom, the study, the kitchen and classroom. It can be standing or wall mounted. The high gloss framing is good looking to suit your home environment perfectly.

8 Languages

The time display and menus can be displayed in one of 8 available languages; English, French, German, Dutch, Italian, Polish, Spanish or Portuguese

Time Mode / Date Mode

The time can be displayed in 12 hours or 24 hours format. The date can be displayed as Month-Date-Year or Date-Month-Year

User Interface Feature

The clock automatically switches between 3 displays; 1. standard mode displaying the Day, Time, Date, Month and Year, 2. Colourful weather forecast, and 3. standard weather forecast.

Alarms

The latest model has; • up to 5 basic alarms (Monday to Friday, Saturday to Sunday, or Monday to Sunday) and •up to 3 Medicine Reminders (one time only, or repeating daily). The user can stop the alarm by pressing any button

Automatic Brightness

The clock has the function to automatically adjust the display depending on the time of day. High, Medium, Low or Auto. In Auto , the display brightens during Day Mode after 7am (50cd/m²) and automatically dims during Night Mode after 7pm (10cd/m²).

Features;

- 8 inch LED Backlight Digital Display Clock
- Day, Time, Date, Month and Year displayed in full wording - non abbreviated. Optional display in analogue clock or coloured background
- Automatic brightness adjustment setting; Day Mode after 7am (50cd/m²) and Night Mode after 7pm (10cd/m²)
- Adjustable volume
- New Model - has up to 5 basic alarms (Monday to Friday, Saturday to Sunday, or Monday to Sunday)
- Up to 3 Medicine Reminders (one time, or repeating daily)
- The user can press any button to stop the alarm
- Display in one of 9 Languages; English, French, German, Dutch, Italian, Polish, Spanish, Portuguese, or Finnish
- 12 hour or 24 hour time format
- Date display format; Date-Month-Year or Month-Date-Year
- Resolution 800 x 600 high definition display
- Display Ratio 4:3
- No glare within 15 feet
- Dimensions 217mm W x 172mm H x 25mm D
- Approximate dimensions of the characters; °Day characters - 25mm H x 12.5-18mm W
- Time numerals - 38mm H x 20mm W
- Date characters - 10mm H x 8mm W
- Weight 410g
- Accuracy +/- 30-60 secs p/mth
- Electricity operated only. Internal battery retains the clock time in the event of a power failure. However, the clock must be plugged in and powered by electricity to operate.
- 5V/1A Power Adaptor supplied

RRP AU \$169



12.1" Digital Calendar Day Clock - Orientation Clock - TTC-DC1201

Orientation Clock for Alzheimer's & Memory Impairment
Non Abbreviated Display of Day, Time, Date, Month and Year
Perfect for people with memory impairment, New Model - Has up to 12 alarms per day



Reduce the workload of carers

The Dementia Day Clock is aimed at helping people living with Dementia by simply and clearly displaying the time of day as morning, afternoon, evening or night. People living with Alzheimer's or Dementia can lose the ability to recognise what the time is and even whether it is day or night. This can cause many challenges for them and for their primary carers. By using this digital calendar orientation day clock it can reduce the challenges.

Help people with dementia live in a daily routine

The aim of the day clock is to help people with dementia to be able to link the time of day to their daily routine. The Day Clocks' clear and simple display of Time, Time of Day and Date, has proven to be helpful for people with dementia or Alzheimer's

Application Scenario

The Day Clock can be used in the living room, the bedroom, the study, the kitchen and classroom. It can be standing or wall mounted. The high gloss framing is good looking to suit your home environment perfectly.

8 Languages

The time display and menus can be displayed in one of 8 available languages; English, French, German, Dutch, Italian, Polish, Spanish or Welsh

Time Mode / Date Mode

The time can be displayed in 12 hours or 24 hours format. The date can be displayed as Month-Date-Year or Date-Month-Year. The Day, Time, Date, Month & Year are all shown, non-abbreviated. Additionally, if Advanced mode is selected in the settings, then the time of day will be displayed i.e. Morning, Afternoon, Evening, Night

Automatic Brightness


The clock automatically adjusts the display depending on the time of day. 200cd/m² auto changed brightness; Day mode: 50cd/m²; Night mode: 10cd/m². You can also set the brightness manually from level 1 to level 5.

Up to 12 Alarms

You can program up to 12 alarms per day. The alarms will repeat daily (or on specific days of the week). The alarms can have one of the following messages on screen at alarm time; Meal, Medicine, Clinic, Shower, Gardening, Cleaning, Appointment and Walk Dog. The user can press any button to stop the alarm.

Features;

- 12.1 inch LED Backlight Digital Display Clock
- Day, Time, Date, Month and Year displayed in full wording - non abbreviated
- If advanced mode selected the time of day is shown i.e. Morning etc.
- Automatic or manual brightness adjustment
- New Model - has up to 12 alarms per day (or on specific days of the week)
- Display in one of 8 Languages; English, French, German, Dutch, Italian, Polish, Spanish or Welsh
- 12 hour or 24 hour time format
- Date display format; Date-Month-Year or Month-Date-Year
- Resolution 1024 x 600 high definition display
- Display Ratio 16:9
- No glare within 15 feet
- Dimensions 299mm W x 207mm H x 27mm D
- Weight 968g
- Electricity operated only. Internal battery retains the clock time in the event of a power failure. However, the clock must be plugged in and powered by electricity to operate.
- 12V/2A Power Adaptor supplied



RRP AU \$179



VISO10 - TTC-VISO10 or TTC-VISO10-WO Large Easy Read Alzheimer's Memory Clock

**Clear Easy to Read - Alzheimer's , Dementia,
Visual and Memory Aid Clock**

*Non Abbreviated Display of Day, Time, Date, Month and Year
Perfect for people with memory impairment*

A digital clock with a non abbreviated display of day, time, date, month and year in written text and numerals. Settings can be locked to avoid tampering. The clock is powered by two AA batteries and has eight language functions. It can be positioned on a table or wall mounted. Suitable for people needing assistance with orientation.

- Time and date display
- Day of week and month display
- Large clear easy to read display
- Non abbreviated display of Day, Date, Month and Year
- 8 languages for day and month display (French, English, German, Italian, Spanish, Dutch, Swedish, Polish)
- 12hr (AM/PM indicator) or 24hr (Military style) display options
- Uses 2 Batteries AA (included)
- Stand on table or hang on the wall
- Fixture hanging hole and table stand on back
- Lock function to disable accidental time change
- Dimensions (mm) : 280 x 190 x 27
- Large Numeric Display (mm) : 190 x 105
- Warranty 2 years from product source

Stylish and functional the new VISO10 Large format Easy Read Alzheimer's Clock is now available from TabTimer Reminders™.

The clear, easy to read VISO10 clock, is a useful tool for people living with memory impairment, caused by Alzheimer's, dementia, stroke, or other memory related impairments - as well as a welcome relief for their family and care givers.

People with Alzheimer's, dementia, stroke, or other forms of memory loss, often lose the ability to recall; what day, time of day, month, or what year it is - causing them stress, anxiety and frustrations. Many large format clocks also abbreviate the day and month, shortening words to three characters (i.e. MON or DEC), subsequently confusing some elderly, or baffling cognitively impaired people.

The VISO10 Clock - now available from TabTimer Reminders - is different! It uses full words, in large bold letters, with no confusing abbreviations. The non abbreviated display of DAY, DATE, TIME, MONTH and YEAR are easy to understand, with large bold numerals up to 190mm high. The battery operated VISO10 also has a choice of 8 different languages, can be desk stand or wall mounted, and can be locked to avoid any unintended time changes. At a total size of 280mm x 190mm the VISO10 is clear and easy to read, whilst remaining sleek and modern, to fit practically into any environment.

The VISO10 Alzheimer's Clock is a perfect gift for home, office, schools, aged care facilities, or the reception area of any health care professional - enabling patients and clients to clearly see the day, date and time - providing a stress free, dementia friendly environment. It may even alleviate carers' frustrations who will no longer have to continually repeat the day and time to their loved ones.

Please Note: This clock does not have any alarms; it is merely a visual and memory aid for the day, time, date, month and year.

The clock also has Atomic Clock functionality which is not available in Australia. But works perfectly fine without it. When used in Europe the clock automatically synchronises daylight saving with the Atomic Clock. As Australia does not have an Atomic Clock transmitter this feature is not available and needs to be adjusted manually for daylight savings changes.

**DVA
FUNDING**

Order now from;

Allianz

RRP AU \$189



**DVA
FUNDING**

Order now from;

Allianz

MemRabel 2 - TTC-MEMRABEL2 Audio Visual Calendar Clock with Reminders

Dementia friendly Audio Visual Calendar Alarm Clock

Perfect for people with Alzheimer's, dementia, memory impairment, cognitive impairment, or everyday forgetfulness.

The NEW MemRabel 2, the latest audio visual reminder alarm clock, for people living with memory impairment, cognitive impairment, or everyday forgetfulness, is now available in Australia from TabTimer Reminders™.

Perfect for people living with Alzheimers, dementia, or other chronic illness or disabilities, the Memrabel 2 is a visually appealing calendar clock with a high quality, bright contrast display, showing the full time, day, time of day, date, month and year, to avoid those repetitive questions about 'what time is it'?

Not only is the MemRabel 2 a functional memory aid clock, but the family or health care worker can easily program up to 20 alarms, with accompanying audio visual and video cues. At alarm time the selected voice reminder with accompanying screen graphic or video is played back on the screen. The alarms can be set as individual one off events, or repeat daily, weekly, monthly or annually, for regular medicines, appointments, birthdays, anniversaries, or any daily living requirement that needs a gentle reminder. The possibilities are endless.

The device comes with over 100 different audio, image and video files pre-installed for typical reminder events, you can create and add your own to the library, or play directly from an external source via USB or SD Card. Reminders can be set to repeat as many times as you think is necessary for each alarm, or repeat at intervals from 5mins to 25mins.

With 4 different calendar display options and four different coloured background options the MemRabel 2 improves visibility for people with poor sight, vision impairment and reduced cognitive ability.

Memrabel-2 has proven successful in helping people live more independently, and has provided families with some peace of mind knowing that in their absence there is a simple assistive technology that really is making a world of difference.

- Large clear easy to read display - time, day and date clock
- Non abbreviated display of Day, Date, Month and Year
- Choice of 4 time display formats
- 12hr (AM/PM indicator)
- Choice of 4 backgrounds and contrast colours
- 20 programmable daily, weekly, monthly, yearly alarms
- Reminders with accompanying screen graphic /photo
- Over 100 pre-loaded voice reminders and graphics
- Choose from pre-loaded voice and picture reminders for medication, eating, drinking, appointments and usual daily chores, or create and install your own
- Infrared remote control included for setup on screen
- Stand on table or hang on the wall
- Powered by 5V mains adaptor with UK>AU adapter (included)
- Long life on-board battery for clock and alarm settings
- Dimensions (mm) : 215 x 175
- Screen size 182 x 140 x 6.4mm, 8" TFT LED Display
- High Resolution 1080/720, Internal file storage 8Gb
- Picture formats JPG, BMP, PNG, GIF
- Audio formats MP3, WMA, OGG, APE, AAC, WAV, FLAC
- Video formats RM, RMVB, AVI, MP4, MOV, VOB, PMP, FLV
- Interfaces Standard USB, Mini USB, 3-in-1 Card Slot SD/MMC/MS
- Warranty 12 months from product source

RRP AU \$299

MemRabel 3i - Touch screen, Wi-fi, APP enabled Orientation Clock TTC-MEMRABEL3i
Dementia friendly touchscreen Wi-fi & APP enabled Audio Visual Orientation/Calendar Reminder Alarm Clock
Perfect for people with Alzheimer's, dementia, memory impairment, cognitive impairment, low vision, or everyday forgetfulness.

MemRabel 3i - touch screen memory prompting aid with wi-fi connection for remote programming

Smart Care from your Smartphone

- ✓ Transfer media files for alarm time playback - once only, Daily, Weekly, Monthly
- ✓ Use personalised videos or photos with text to create memory prompts or routine planning
- ✓ Easy to use music and video recording player
- ✓ Create and send media file memory prompting alarms with your Smartphone

The MemRabel 3 is a high-performance prompting calendar clock. MemRabel 3 uses touch-screen technology to provide a simple method of creating personalised memory prompts (unlimited alarm reminders) to benefit people with cognitive impairment. Videos and photos are easily transferred using an SD card or USB Flash Drive to use as a memory prompt. You can split the screen and add personalised text to accompany video or photo playback at alarm time.

User interaction level with MemRabel 3 can be selected:

- No interaction alarm: memory prompts are played automatically at selected times. Screen is passcode locked to prevent tamper.
- Video playback: Upload video files and assign to touch-screen playback thumbnail tiles.
- Music player: Upload a favourite music playlist for simple one-touch play/pause/resume/stop. All other settings are locked out.

MemRabel 3 allows carer instruction videos to be created for playback. The carer can play a family created video to learn about new medications or care preferences eliminating the need for whiteboard or written messages.

MemRabel 3 can be connected to home or business Broadband to allow memory prompt files to be created and sent over the internet to set up reminder alarms with associated media files. The Smartphone APP is available for iOS and Android and can be securely used by multiple family/carers members.

Key Features

- Timed playback of personalised video, personalised photos with text, or slideshows
- Touch screen alarm time memory prompt set up
- Plan a daily routine with video instruction
- Manage medication reminders
- Select alarm repeat for single, daily, weekly or monthly
- Create and control timed playback of family or memory prompting slideshows
- Choice of clock face and colour
- Daily repeat medication course setting (i.e. for 8-10 days)
- Easy music player upload a favourite playlist of music
- Assign your videos to a touch tile for easy playback
- Create carer instruction video for one-touch playback
- Connect to internet for free software updates (optional)
- User friendly Smartphone APP
- Simplified alarm memory prompt set up wizard
- Create and send messages or reminders from anywhere
- Media file copy wizard from Flash Drive
- Memory prompt alarm files arranged for easy view and edit
- Night mode display dimmer, to select screen dimming time
- Touch screen music player - play/pause/resume/shuffle/stop
- Touch screen memory prompt media file playback stop
- Hearing impaired visual alarm notification
- Alarm countdown alternates screen from black to white with numbers descending from 10 to 1 preceding alarm playback
- Choice of user interaction level with screen passcode lock
- Automatic seasonal time adjustment
- Use a single phone with multiple devices (Android only)
- Use multiple phones with a single device (for family)
- Over 100 stock reminder files pre-installed prompt creation
- Automatic free software updates when connected to home or business Broadband
- No subscription charges, contracts or agreements to sign (Over the internet service is free for the first 24 Months. Thereafter buy 10 message credits for \$10)
- See usage policy, privacy and personal data notices on our website

RRP AU \$399



VISO 50 - Analogue clock and digital orientation display –TTC-VISO50

Alzheimer's Visual & Memory Clock

Perfect for memory impairment



Clear Easy to Read Analogue Clock with Non-Abbreviated Day and Month, Time and Date Display, Day of week and month display

Analogue clock (200×200 mm) : hours / minutes / seconds hand

Digital display clock (200×90 mm), 4 display formats (non-abbreviated words) :

1. Day of the week / Part of the day - i.e. Friday Afternoon
2. Month / Day of the month (digit) / Year - i.e. January 14 2019
3. Day of the week / Hour - i.e. Friday 10:00
4. Adjustable 3 display formats switch from every 5 to 60 seconds

- 12/24H display for the digital display
- 5 languages available : French, English, German, Italian, and Spanish
- Wall or table mount
- Dimensions : 260 x 380 x 45 mm
- Uses 3x AA batteries (not included)
- Fixture hanging hole for wall or built in stand

RRP AU \$210

Motion Activated Memo Reminder-TT-MAMR

Personalised recording prompt for memory impaired

Remind your loved one when they activate the sensor



When people are a little forgetful, this movement activated voice memo player can help by using family recorded messages or instructions. You simply record a voice message then switch the alarm on. When a person is detected by the photoelectric beam the memo is played back.

Use to remind an elderly person not to let strangers into the home or to check a callers identity, to warn children not to go outside, or as a welcome message in reception to greet guests.

- Powered by 3 x AA batteries
- 3-minute delay between voice announcements (prevents repeats if a person stands in the beam)
- Detection range 2-3 metres
- On/off/record switch
- Swivel wall bracket for optimum detection angle
- Supplied with all fixings and batteries.

RRP AU \$99

Assistive Listening Devices

Artone 3 Max – ART-3MAX

World's Strongest Streamer Bluetooth loopset/neckloop works even for severe hearing loss and cochlear implants.

Artone-3 MAX streamer is the most powerful bluetooth neckloop on the market that works even for severe hearing loss and cochlear implants.

This version of the popular Artone bluetooth loopset/neckloop will create a link from your hearing aid to your Bluetooth phone. When purchasing the Artone-3 product you upgrade your hearing device to 'Bluetooth hearing aid' level. Hearing aid/implant with Bluetooth function allows the user to communicate over cell phone, using the hearing aid as a handsfree.

Further more you can listen to the TV while using the hearing aid/implant as wireless headphones. Artone-3 streamer also allow music streaming from iPhone, iPods, iPads straight into hearing device. Artone 3 users can enjoy: podcasts, videos, music transmitted straight to the hearing aid/implant.

SHORT DESCRIPTION

ARTONE 3 streamer creates a link from a hearing aid/implant (must be equipped with T-Coil) and any Bluetooth-enabled device such as mobile phone, computer, music player etc. Thereby, ARTONE 3 enables crystal clear sound transmission for hearing aid users during mobile phone conversation or music streaming.

WHAT ARTONE 3 DOES

ARTONE 3 upgrades a T-Coil equipped hearing device in to a Bluetooth hearing device.

Using the ARTONE 3, the hearing aid users can communicate over mobile phone using their hearing aid as handsfree or listen to music using their hearing aid as wireless headphones.

Crystal clear amplified sound of voice or music is transmitted from the audio source to the hearing aid through the ARTONE 3 without any interference.

HOW ARTONE 3 WORKS

1. USING THE ARTONE 3 WITH A MOBILE PHONE:

In a time of a call the signal from the mobile phone is wirelessly sent to the ARTONE 3 Bluetooth LoopSet where the signal is processed, cleared and amplified. Then the signal is sent to the user's hearing aid or implant.

The hearing aid processes the signal once again and adjusts the signal to the user's personal level. When having a mobile phone conversation, the hearing aid user's voice is captured by the ARTONE3 Bluetooth® LoopSet microphone, it is processed and cleared by cutting the unnecessary background noises, it is amplified and then the signal is sent to the user's handset.

The product is equipped with noise cancellation system that allows the hearing aid user to communicate even in a noisy environment.

2. USING THE ARTONE 3 WITH BLUETOOTH-ENABLED AUDIO SOURCES

When paired with any audio device equipped with Bluetooth, ARTONE 3 enables the transmission of crystal clear amplified sound of music, film, Skype conversation etc. without any interference from the audio device into the hearing aid.

To enjoy crystal clear amplified sound streaming from non-Bluetooth audio sources such as television, Ipod, MP3 player etc. or for amplification of conversation in noisy environment, an additional transmitter must be used with the ARTONE 3 – the ARTONE TVB transmitter or the ARTONE TVB+ Wireless Microphone and Transmitter.

KEY FEATURES

- Supports Bluetooth A2DP, AVRCP, and HS/HF profiles
- Wirelessly enjoy stereo music from all A2DP-enabled devices, such as iphone cell phones, PDA, PC etc
- Features Software & CVC echo/noise cancellation for reducing noise
- DUAL MICROPHONE, ECHO & NOISE CANCELLATION
- Remote controls: play, skip and select your music.
- Totally wireless.
- Clear amplified communication.
- Compatible with any phone with Bluetooth® technology.
- Compatible with any hearing aid equipped with T-coil option Rechargeable Long lasting Battery
- Voice dial options.
- compatible with Artone TVB – for SKYPE, TV, Radio, Ipod...)
- Light and elegant style.

people are a little forgetful, this movement activated voice memo player can help by using family recorded messages or instructions. You simply record a voice message then switch the alarm on. When a person is detected by the photoelectric beam the memo is played back.

Use to remind an elderly person not to let strangers into the home or to check a callers identity, to warn children not to go outside, or as a welcome message in reception to greet guests.

- Powered by 3 x AA batteries
- 3-minute delay between voice announcements (prevents repeats if a person stands in the beam)
- Detection range 2-3 metres
- On/off/record switch
- Swivel wall bracket for optimum detection angle
- Supplied with all fixings and batteries.

RRP AU \$210 \$199

VibraLITE 12 Mini **VibraLITE 12 Mini - Vibrating Reminder Watch**

Up to 12 Daily Alarm Settings

Smaller Case & Watch Band - For Small Wrists of All Ages



Model TTW-VM-LBK / LPL

AU \$109

(leather buckle style band)

Black Fits Wrists: 11.3cm to 18.25cm

Purple Fit Wrists: 11.8cm to 18cm



Model TTW-VM-SWH / SMC

AU \$109

(Silicon buckle style band)

Fits Wrist Sizes: 12.5cm to 18.4cm



Model TTW-VM-VBL, TTW-VM-VOR

Model TTW-VM-VPN

AU \$109

(Velcro hook and loop style band)

Fits Wrist Sizes: 13.75cm to 18cm



Model TTW-VM-FLEX AU \$129

(Stainless steel flexi stretch band)

Fits any wrist add/remove links to fit

- Smaller Case (30.5mm x 38.5mm) (12.3mm Thick)
- Easy to read time digits (5mm H x 3mm W)
- Attractive Styling - 5 New Models to Choose From
- Water Resistant 20 meters (5 ATM)
- Accuracy per month: 30 seconds (+/-)
- Automatic Leap Year Handling
- One Year Limited Warranty from source (excludes battery)

New Lock-Out Feature:

Press the "MODE" & "SET" buttons at the same time while in "TIME" mode to lock/unlock the watch, so time/countdown timer/alarm settings cannot be changed. The EL (electro-luminescent) light will still function, as will the ability to manually stop alerts (to conserve battery).

ALL VibraLITE 12 Mini Watches have the features below

Vibration, Sound or Both

- Securely set in a separate "OPTION MODE"
- 20 second alarm duration with auto off or push any button

12 Daily Alarm Settings - Preset for variable reminders

- Examples: 7:00 am - wake up
- 7:30 am - first daily medicine
- 8:20 am - time to leave for work
- 10:30 am - daily medicine due etc

1 x Countdown Timer with Repeat Option

- For ongoing preset interval timing
- Alarm sounds when countdown duration reached
- Repeat option to repeat the countdown
- Counts down in hours, minutes and seconds- up to 23 hours and 59 minutes and 59 seconds

EL Back Light

- Stays on while any button is being pushed making it easier to set in the dark

Text Prompts

- Assist in setting and confirmation of settings

Large Digits Makes for easy reading

12 Hour Clock Uses am/pm clock

RRP AU \$109-129

**DVA
FUNDING**

Order now from;

Allianz

e-pill CADEX VibraPlus Vibrating Alarm Watch



Black Leather Band
TTW-CAD-VP8
RRP \$179

Backlit Display

Easy-to-read display shows time, mode, and alarms set.



Lock Feature

Buttons can be locked to prevent accidental re-programming.

Leather Strap

This watch uses a standard 18 mm replaceable strap.

Up to 8 Alarms

When alarms are set they appear on the display. They can be set to sound, vibrate, or both.



Stainless Steel Flexi Stretch Band
TTW-CAD-VP8-FLEX
RRP \$199



Lanyard Kit Pendant Conversion

TTW-CAD-PEND
RRP AU \$ 19

e-pill CADEX VibraPlus 8 Alarm Vibrating Watch

Up to 8 Daily Alarm Settings with Vibration
STRONG Vibration Alarm, Sound Alarm, or Both

The e-pill CADEX VibraPlus is a nicely designed medication reminder alarm watch with strong vibration or sound. It looks like a normal watch; it shows real time, up to eight daily alarms that auto-reset each day at midnight (repeats again tomorrow at the same times), has a back-light, stop watch and unique Count-Down Timer with "auto start" (interval timer automatically starts and stops at the end of each day).

Watch Vibrates! - time to take your pills, get an injection, or remember that important activity!

- Multi-Alarm (up to eight alarms per day)
- Discreet & Easy to Use Medication Pill Reminder Alarm Watch
- Black/Silver polycarbonate case (40mm X 50mm)
- Display Time Digits (6.35mm H X 3mm W)
- Stainless Steel Back
- Replaceable black leather standard 18mm strap
- Vibration Alarm - discreet alarm helps hearing impaired patients
- STRONG Vibration Alarm, Sound Alarm, or Both
- Count Down Interval Timer - long interval from 15 minutes to 24 hours with unique Auto Start
- Key Lock feature (prevents accidental changes)
- Back Light makes it easy to read the time at night
- Stopwatch (Chronograph)
- Battery: Replaceable CR2032 3 Volt, Battery life 4-6 months.
- 1x Spare Battery and screwdriver included to facilitate change
- Easy to Read Instruction Manual
- Water Resistant
- One Year Warranty from product source (excludes battery)

This watch is for anybody who does not want the world around him or her to know they are on medication.

"This watch has allowed me to be prompt with my eight times a day pill taking. This compliance helps me stay "steady" all day. In the past I either forgot or was late or couldn't remember if I had taken my meds. I am really pleased with the watch that vibrates on my wrist and doesn't announce to the rest of the world that it is "time to take your pill." P.V. NJ

Originally developed for children with diabetes - e-pill CADEX VibraPlus watch vibrates (strong vibration) - time to go to the school nurse & get an insulin injection! Not only for children - this watch style is for adults and teenagers! Great looking watch for individuals of all ages.

The e-pill CADEX VibraPlus Alarm Watch works as a discreet reminder in meetings, school, church etc. Size is suitable for men, women and children.

Never forget to take your medications again.
Avoid medication errors. Feel better.

RRP AU \$179 – \$199

Watch Minder Watches

WatchMinder 3 - vibrating watch reminder system Up to 30 Daily Alarm Settings with Vibration STRONG Vibration Alarm

Choice of Black, White, Blue, Red or Pink Silicone Band with Stainless Steel Buckle. Fits wrist sizes 12.5cm to 20.5cm

Simple. Programmable. Discreet.

The vibrating wristwatch that will privately cue you through your day.

Choose a pre-programmed message or create your own personalised message. Choose a fixed time or interval. The WatchMinder will discreetly cue you through your day.



WM3-BK



WM3-WH



WM3-BLU



WM3-RED WM3-PK



MS-4007 \$ 19



MP3310 \$ 29

- Black polycarbonate with fibre case (38mm diameter)
- Stainless Steel Back
- Real Time, Calendar Clock & Date
- Vibration Alert Only, No audible alarm
- **Training Mode** - one message at intervals throughout the day
- Training mode alert, lasts 10 seconds, repeats in 20 seconds, push any button except light to stop. e.g. RELAX every 15 min
- **Reminder Mode** - up to 30 reminders per day, repeat daily
- Reminder mode alert with snooze, alert lasts 10 seconds, repeats in 20 seconds. Push any button except light to stop the vibration
- **Reminder Snooze Feature** automatically repeats the first cycle in 2 minutes
- Create your own 12 Personalised Messages up to 8 characters
- Display Time Digits (6mm H X 2mm W)
- 12 hour clock display, AM/PM
- Weighs 33 grams
- Waterproof
- One Year Warranty from product source (excludes battery)
- Charged via provided USB cable, plug it into any USB port to charge
- Need a USB port - purchase optional MS-4007 or MP3310
- Fits wrist sizes approx. 12cm to 20.5cm

The WatchMinder is a simple wristwatch that can easily be programmed to set up discreet vibrating reminders throughout one's day. Invented by a child psychologist, the WatchMinder was originally designed to aid children with Attention Deficit Disorder (ADHD). It is now used by people of all ages, varying chronic illness and disabilities, and others with special needs in staying focused, managing their time, remembering meds or tasks and modifying their thoughts and behaviours.

- Attention Deficit Hyperactivity Disorder (AD/HD)
- Brain Injury & Stroke
- Self Monitoring
- Learning Disabilities (LD)
- Medication Compliance
- Independence Training & Vocational Rehabilitation
- Elderly
- Deaf and Hard of Hearing
- Research, Data Collection & Clinical Trial Studies
- Chronic Illness
- Behaviour Modification
- Assistive Device

There are 15 preset training mode messages; Bathroom, Be Positive, Breathe, Cough, Give Positive, I can Do It, Ignore, Pat ATTN, Posture, Pray, Relax, Slow Down, Smile, Stop ,Stretch

There are 50 preset reminder mode messages; Backpack, Bank/Bills, Bathe, Bathroom, Bed Time, Books, Call/Text, Chores, Clean, Cook, Diabetic Test, Doctor, Dress, Eat, Exercise, Get Kids, Go Home, Go Lesson, Go Meeting, Go School, Go Sport, Go Work, HW/Study, Mail/Email, Meds 1, Meds 2, Meds 3, Meds 4, Meds 5, Meds 6, Meds 7, Meds 8, Pet Care, Planner, Play Time, Pray, Read/Write, Recharge, Rest, Snack, Store/Gas, Stretch, Talk Teacher, Teeth, Trash, Turn in HW, TV/DVR/DVD, Water, Write, Write HW

RRP AU \$179

Low Vision & Braille Talking Watches



TTW-LVTW-WBK \$59

Low Vision Talking Watch - TTW-LVTW for low vision or vision impaired (SOUND ALARM & TALKING VOICE ONLY, DOES NOT VIBRATE)

Speaking watch

- Alloy made watch case with nickel free processing
- Watch face diameter 40mm
- White watch face with large bold black font
- Watch hands with contrast colour
- Japanese quartz movement
- Genuine leather band (Black, White, Brown or Yellow)
- Water resistant watch bottom case
- Depth of watch 12MM. total length of watch 240mm
- User manual and gift box packaging
- Weight 52g



**TTW-LVTW-WHT
\$59**



**TTW-LVTW-BRW
\$59**

RRP AU \$59 - \$75



**TTW-LVTW-YELL
\$59**



**TTW-LVTW-GLD
\$75**

Braille Low Vision Talking Watch - TTW-BRTW-WBK for blind people, low vision or vision impaired (SOUND ALARM & TALKING VOICE ONLY, DOES NOT VIBRATE)

Speaking watch

- Alloy made watch case with nickel free processing
- Watch face diameter 40mm
- White watch face with large bold black font
- Three dots indicate 12 o'clock; two dots indicate hour 3, 6 and 9; one dot indicates hour in between
- Lid opens at 6 o'clock position
- Professional dot design for excellent touch feeling
- Swiss quartz movement function
- Black Genuine leather band
- Water resistant watch bottom case
- Depth of watch 12MM. total length of watch 240mm
- User manual and gift box packaging
- Weight 64g



TTW-BRTW-WBK

**PLEASE NOTE:
has leather watch strap,
not stainless steel**

RRP AU \$69

Vibrating Clocks



Sonic Alert Sonic Shaker™ - SBP100 - Sonic Boom - Vibrating Travel Alarm Clock

Small, light-weight, self-contained and battery-operated, this is a vibrating travel alarm clock that you can count on
The SBP100 Sonic Boom Alarm Clock is small enough to fit in the palm of your hand, but packs a pretty big wake-up punch!

- 90 db Extra-loud alarm
- Sonic Shaker™ bed vibrating unit
- Easy to set function keys, Easy to read .75" display
- Continuous four minute Snooze setting
- New green power technology now uses 60% less energy
- Standard size batteries included (2 AA; 1 AAA)
- Travel case with pillow clasp included
- One-year Sonic Alert® blast-proof warranty

RRP AU \$88

ILY Clock with Motivator Bedshaker - TTC-ILY

WAKE UP WITHOUT DISTURBING OTHERS

Loud Quartz Analogue Alarm Clock (Plain face) - up to 96dB with Motivator Bedshaker / Pillow shaker



FEATURES:

- Plain Face, Easy to use analogue style
- Variable Tone (select low to high) ,
- Variable Volume (up to very loud 96 dB)
- Snooze Function, Button Activated Night Light
- Vibrating bed shaker included
- Battery Operated (Requires 4x AA Batteries, not included)
- Optional power supply extra
- One Year Warranty

USES:

- Wake to Sound, Vibration or Both
- Great for Travel; Runs on Batteries . . . or with optional AC Adapter
- Great for waking up . . . without disturbing others!



Suggested for:

- Deaf & Hard of Hearing
- Heavy Sleepers, Shift Workers , Early Risers
- Students, ADHD , Elderly , Travellers

RRP AU \$99

TIME TIMER®

Visual Timers

For
**Home
Work
Education
or
Special Needs**



Time Timer PLUS® 5 Minute – TT-TTP05-WHT

Up to 5 minute visual countdown timer

The best tool for keeping short routines at home, school or the office on track

The Time Timer PLUS 5 Minute offers a greater visual impact for shorter duration activities such as brushing teeth, holding a speed round meeting, breathing exercises and managing a brainstorming session. The Time Timer PLUS 5 Minute is the best tool for keeping short routines at home, school or the office on track.

Benefits

- Teach children the concept of elapsed time
- Perfect for short-duration activities
- Hold a speed round meeting
- Reduce conflict in family and household routines
- Manage the stress of transitions by showing “how much longer”
- Empower teachers to teach instead of constantly monitoring behaviour
- Transcends language and culture barriers

Features

- Volume control dial for those times you need a louder tone and for those times when you prefer silence
- Easy-to-read numbers, Lightweight anti-static disk
- No loud ticking, Size: 13.75cm x 17.5cm
- One AA battery required for operation (included)
- Made in China

RRP AU \$99



Time Timer PLUS® 60 Minute– TT-TTP60-CHA

Up to 60 minute visual countdown timer

Fast becoming the most popular Time-Timer for teachers, students, parents and professionals, with on-the-go handle makes it easy to take anywhere

Fast becoming the most popular Time-Timer, the Time Timer PLUS® 60 Minute packs in all the features that have made Time-Timers the first choice of teachers, students, parents and professionals into a new design that's as sturdy as it is useful.

Ideal for active kids or a day full of meetings, the Time Timer PLUS 60 Minute features a rugged case and durable clear lens to protect its patented red disk. And the new on-the-go handle makes it easy to take the timer anywhere you go!

Benefits

- Teaching children the concept of elapsed time
- Making homework and practice time more productive
- Reducing conflict in family and household routines
- Managing the stress of transitions by showing “how much longer”
- Empowering teachers to teach instead of constantly monitoring behaviour
- Increasing performance on tests and in the workplace
- Keeping meetings on time and moving forward
- Transcending language and culture barriers

Features

- Volume control dial for those times you need a louder tone and for those times when you prefer silence
- Easy-to-read numbers, Lightweight anti-static disk
- No loud ticking, Size: 13.75cm x 17.5cm
- One AA battery required for operation (included)
- Made in China

RRP AU \$99

Featured Award



This item is recognized by Able Play - a rating system that focuses on the developmental areas of physical, sensory, communicative, cognitive, and social/emotional abilities of a child.

Time Timer® MOD + Dry Erase Board – TT-TTMB60

Up to 60 minute visual countdown timer with marker board

Style and durability in a portable size that's easy to use anywhere



The Time Timer MOD + Dry Erase White Board includes both a Time Timer MOD in Charcoal and the Time Timer Dry Erase Board.

The Time Timer MOD is a 60-minute visual timer that has the style and durability of the Time Timer PLUS in a size that's portable and easy to use anywhere—from your office desk to the kitchen counter. The Time Timer Dry Erase Board is an activity and message centre for office, home and school. It is perfect for time-sensitive messages and keeps schedules on time.

Dry Erase markers not included.

Benefits

- Stay on track by showing exactly what needs to be completed in the time available
- Ensure each step in a routine is followed - nothing is missed!
- All important topics for a meeting are covered in the time allowed

Features

- Convenient caddy to hold markers and supplies
- Size: 20cm x 20cm, Clock : 9cm x 9cm
- One AA battery required for operation (included)
- Made in China & USA

RRP AU \$120

Time Timer 90-minute TWIST® - TT-TTT90

Up to 90 minute visual countdown timer

The latest time management tool from Time Timer. Don't settle for imitations

The 90-minute TWIST® is an easy-to-use time management tool great for kids and adults. To use, simply turn the outer ring until you get to your desired time, then push the play button. A digital display shows at a glance how much time remains. A magnetic back makes the TWIST® a snap to display. Easily attach it to most metallic surfaces, from refrigerators and file cabinets, to the whiteboard in your classroom.

The 90-minute Time Timer TWIST® is a quiet and versatile timer that helps make all kinds of everyday tasks easier and more productive:

- At home use the TWIST® for everything from keeping kids' morning and bedtime routines on task, to timing homework, instrument practice and screen use
- In the classroom, it's great for everything from timing tests and turn-taking exercises, to timing group activities and helping kids with special needs stay on task
- In the office, use it to keep track of presentations, team-building activities, and more

Benefits

- All it takes is a quick twist of the red outer ring to set this timer for up to 90 minutes, as required, then press the play button
- A magnetic back makes it easy to put visual time management right where you need it. Perfect for the office, classroom or kitchen
- Keep the Time Timer TWIST on the whiteboard in your classroom so that it is ready to use at any moment

Features

- Play / Pause button for quick starts and easy pause functionality
- Digital tick marks appear when timer is actively running
- Duration of timer appears in the centre
- Quiet operation, no loud ticking
- Digital alert sounds at the end of the timer
- Magnetic back to allow for easy placement
- Size : 8cm diameter x 2cm H
- One AAA battery required for operation (included)
- Made in China



RRP AU \$66



Featured Award



This item is recognized by Able Play - a rating system that focuses on the developmental areas of physical, sensory, communicative, cognitive, and social/emotional abilities of a child.

Time Timer Watch PLUS® Small – TT-TTWSML-CHA

Watch with Visual countdown timer and alarm

The perfect solution for on-the-go, convenient time management

Watch your time be even more productive with the Watch PLUS®. An easy way to keep every daily routine on track, the Watch PLUS is the perfect solution for on-the-go, convenient time management. The Watch PLUS can help anyone sail through their day, from timing exercise routine intervals to your big presentation. The sporty design makes this wrist-watch a great accessory for students, adults, athletes, teachers and other professionals.

3 Powerful Modes!

- Time Timer: Set an original 60-minute or custom timer up to 99 hours with optional repeating alert
- Clock: Analog and digital displays in 12- or 24-hour format
- Alarm: Set a time-of-day alarm

Benefits

- Silent operation when switching between modes makes this watch perfect for standardized college entrance testing
- Optional vibrating alert makes this watch perfect for meetings or anytime you want discrete time management
- Easy way to time exercise intervals to keep workouts on track

Features

- Vibrating Alert, Repeatable Time Segments
- Backlight, Water Resistant (30M)
- Soft silicone watch band fits wrist sizes between 11.5cm – 18cm
- Watch case body is 3.75cm W x 3.5cm H, LCD screen 2cm x 2cm
- The battery (Type CR2032) is fitted at the time of manufacture
- Made in China

RRP AU \$199

Time Timer Watch PLUS® Large – TT-TTWLRG-CHA

Watch with Visual countdown timer and alarm

The perfect solution for on-the-go, convenient time management

Watch your time be even more productive with the Watch PLUS®. An easy way to keep every daily routine on track, the Watch PLUS is the perfect solution for on-the-go, convenient time management. The Watch PLUS can help anyone sail through their day, from timing exercise routine intervals to your big presentation. The sporty design makes this wrist-watch a great accessory for students, adults, athletes, teachers and other professionals.

3 Powerful Modes!

- Time Timer: Set an original 60-minute or custom timer up to 99 hours with optional repeating alert
- Clock: Analog and digital displays in 12- or 24-hour format
- Alarm: Set a time-of-day alarm

Benefits

- Silent operation with vibrate alert option makes this watch perfect for standardized college entrance testing
- Vibrating alert makes this watch perfect for meetings or anytime you want discrete time management
- Easy way to time exercise intervals to keep workouts on track

Features

- Vibrating Alert, Repeatable Time Segments
- Backlight, Water Resistant (30M)
- Soft silicone watch band fits wrist sizes between 14cm – 21cm
- Watch case body is 4cm x 4cm, LCD screen 2.5cm x 2.5cm
- The battery (Type CR2032) is fitted at the time of manufacture
- Made in China

RRP AU \$199



Voice and sound activated

Voice and Sound Activated Pager – Amplicall 30/50 – TT-VSAP

Sound monitor alerter for deaf and hearing impaired
Baby Monitor, Door Bell, Phone, or Smoke Alarm Monitor for deaf and hearing impaired people. Carer call system for disabled and incapacitated. Use your voice, or sound, to create an alarm.



For deaf & hearing impaired

DEAF BABY MONITOR

DOOR BELL MONITOR

PHONE ALERT

SMOKE ALARM MONITOR

For disabled and incapacitated

CARER HELP CALL ALERTER



AC30 Sound Monitor

- The AC30 is a sound detecting monitor that sends a wireless signal up to 75 meters to the AC50 vibrating beeping pager receiver, when a sound is present.
- Battery operated but can use optional power adaptor.
- Sound sensitivity adjustment and delay control.
- The sensitivity control can be adjusted so that very quiet sounds will cause an alarm, or so that louder sounds are required to cause an alarm.
- The delay control adjusts so that an alarm will be triggered immediately a sound is detected, or can be varied to require sounds lasting up to 30 seconds before an alarm is activated.
- Emergency call push button use also possible
- Jack socket for connection to an external alarm system
- Operates on batteries included 3xAAA

AC50 Vibrating Beeper Pager Receiver

- The AC 50 is a portable lightweight pager powered by included 3 x AAA batteries. The pager has a vibration alert or tone beeper (selectable by switch) and a tone beeper volume switch.
- The pager can indicate by flashing LED, beep and vibration, triggers from up to four paired transmitters, enabling multi-patient or multi-product monitoring. Example: Baby monitor, door bell, personal help call pendant, wander detecting device.

- Easy pairing system
- Pair with up to 4 different devices
- Volume control (off, low, high) up to 80dB
- Vibrating switch on/off
- Low battery indicator
- Test button
- Sensitivity adjustment
- Possibility to connect a plug adaptor power source
- 3.5mm Jack socket for connection to an external alarm

- The AC30 and AC50 together provide quality alarm receivers for home care alerts to carers. The possibilities are endless and are only limited by your imagination.

RRP AU \$169

MedPage Vibrating pager with splash proof Transmitter Bundle - MPPL-ERTXSET

RRP AU \$199



MedPage Vibrating pager with splash proof Transmitter Bundle - MPPL-ERTXSET AU \$199

Medpage Beeping Tone & Vibrating Alert Pager/Receiver

The new MedPage model MPPL carer alert pager offers a versatile range of options for care solutions around the home. Wander alerts, help call carer alerts, pressure mat alarms, door and window alarms, can all be sent to the pager to raise an alarm. With a choice of alarms (vibration only or vibration with sound alarm) the pager is suited to all ages and sensory ability for deaf and hearing impaired or people with disabilities, risk of falls or need to call their carer within the home setting.

The MPPL pager is the key alarm, or alerting device, signalled by the various sensor transmitters. When warnings are needed to alert a carer or partner to daily living scenarios the MPPL pager provides a reliable, economic solution. You can use various transmitters to activate the vibrating pager. Once activated the pager will continue to vibrate and beep (if selected) until the power/reset button is switched off/on.



Medpage Beeping Tone & Vibrating Alert Pager/Receiver - MPPL AU\$139

- Transmitter code learn-in button
- Stylish durable design, Rugged belt clip
- Powerful vibration, Tone beepers with volume selection
- Frequency 433.92mhz
- Operating range 120 metres (open field)
- Dimensions: 80mm x 55mm x 22mm
- Weight including batteries 120g
- Fixing for lanyard, Low battery warning LED
- 3 channel alert (identify different transmitter types)
- Self test at switch on., Powered by 2 x AA batteries
- Latching alarm - must be reset to cancel alarm
- Alerts: vibrate only, vibrate+low tone, vibrate+high tone
- Use up to 10 transmitters with single pager
- Use any number of pagers with single transmitter

RRP AU \$139



MedPage Waterproof Transmitter Help Call pendant - MPP-ERTX AU \$99

MedPage splash proof Transmitter Help Call pendant

MedPage MPP-ERTX (splash proof Transmitter) is a special transmitter housed in a splash proof case with an easy-press switch. The splash proof transmitter is supplied with a lanyard and a hook and loop fastening wrist strap.

The transmitter button can be operated by people with limited dexterity, people with Arthritis or disabilities etc.

Call a carer via wireless technology using the beeping and vibrating alert pager. Help call alerter pager.

Designed for use with the MedPage Beeping Tone & Vibrating Alert Pager [MPPL-PAGER] (purchased separately)

Powered by 1 x A23 battery.
Transmitter range 120M.
Dimensions: H 75mm 50mm D 20mm

RRP AU \$99

TabTimer Vib-Ra-Bell – wireless door bell & vibrating shaker receiver bundle - TT-VB2+1

Designed for the blind, hearing impaired and senior care. Great for noisy environments where the doorbell can't be heard.

Bundle includes; Vib-Ra-Bell – wireless door bell - TT-VB1 and Vib-Ra-Bell - shaker receiver - TT-VB2. The shaker receiver will be paired with the doorbell on arrival.

52 ringtones to choose from and 4 steps volume control (25~110dB), 250~300m working range. Approved with CE, FCC, RoHS certification. Easy to install, no electrical connection or wiring required, since it operates by radio frequency.

When paired with the Vib-Ra-Bell wireless doorbell the Shaker Receiver will flash, vibrate and can emit an alarm sound (25-110 dB) leaving you free to roam around your home without missing that important visitor. Perfect for those that cannot hear the door bell and may be deaf or hard of hearing.



Battery-operated and portable, Mounting screws and adhesive tape supplied, receiver can be worn on waist belt with supplied belt clip, 2x AAA & 1 x LR23 batteries included, Built-in mini vibrating motor with LED indicator, 52 high-quality melodies, 4 adjustable volumes, 250~300m working range, Learning code, can work with different sensors, Power on/off button, Music memory function for each kind of sensor. When the signal is received from the transmitter, the receiver will vibrate, play music and the LED indicator will flash at the same time.

RRP AU \$169

TabTimer Vib-Ra-Bell – wireless door bell TT-VB1

Requires receiving device Vibrating Ringing Indicator - Vib-Ra-Bell Shaker Receiver - TT-VB2. The doorbell can be paired with the shaker receiver very easily.

When paired with the Vib-Ra-Bell Shaker Receiver the alerter will flash, vibrate and can emit an alarm sound (25-110 dB) leaving you free to roam around your home without missing that important visitor. Perfect for those that cannot hear the door bell and may be deaf or hard of hearing.

Battery-operated and portable, Mounting screws and adhesive tape supplied, 1x LR23 battery included, 250~300m range.



RRP AU \$89

TabTimer Vib-Ra-Bell – Shaker Receiver TT-VB2

Requires transmitting device Vib-Ra-Bell – wireless door bell TT-VB1. The shaker receiver can be paired with the doorbell very easily.

When paired with the Vib-Ra-Bell wireless doorbell the Shaker Receiver will flash, vibrate and can emit an alarm sound (25-110 dB) leaving you free to roam around your home without missing that important visitor. Perfect for those that cannot hear the door bell and may be deaf or hard of hearing.

Battery-operated and portable, Can be worn on waist belt with supplied belt clip, 2x AAA batteries included, Built-in mini vibrating motor with LED indicator, 52 high-quality melodies, 4 adjustable volumes, 250~300m working range, Learning code, can work with different sensors, Power on/off button, Music memory function for each kind of sensor. When the signal is received from the transmitter, the receiver will vibrate, play music and the LED indicator will flash at the same time.



RRP AU \$89

Amplcall 2 – emergency alerter / carer alert TT-AC2

Requires receiving device Vibrating Ringing Indicator-AMPLICALL50.

The Amplcall 2 has a long pull cord that can be used when a person has fallen down and needs assistance. It also can be placed in a bathroom so that a person can summon help whilst in the bath or the shower.

The Amplcall receivers all emit extremely loud sounds and have bright flashing lights to attract speedy attention to the emergency.

Product Features:

- White illuminated doorbell with pull cord
- Easy pairing between with the Amplcall 20 & 50
- The Amplcall 2 is perfect for using in bathrooms, kitchens and any place where there is a likelihood that someone will need emergency assistance after falling

RRP AU \$119



With Alert Signaller Amplcall 50 AU \$169

Amplcall 30 – Baby Monitor/ Sound Monitor TT-AC30

The Amplcall 30 is a wireless baby alarm/sound activated transmitter part of the amplcall 30/50 alerting system.

Please note that the Amplcall 50 must be used in conjunction with this transmitter.

RRP AU \$119



Amplcall 16 – Wireless Door Bell Amplifier TT-AC16

Never miss a door bell again with the AmpliCALL16!

When your door bell rings, this doorbell amplifier indicator will ring very loudly and flash to let you know there is someone at the door.

Product Features:

- Extra Loud Ringer (up to 95dB), Adjustable volume
- Camera Flash Ring Indicator
- Battery operated 4x AA supplied, Low battery indicator
- Test button on main unit
- Long wireless range from the doorbell to the main unit(up to 30m)
- Wall and Desk Mountable
- Optional power adaptor MP3310 purchase separately

RRP AU \$169





Amplicall 20 Wireless Doorbell & Telephone Ring Indicator TT-AC20

Never miss a door bell or a telephone call again!

When your door bell/telephone rings, this doorbell/telephone ring indicator will ring very loudly and flash to let you know there is someone at the door/ or that there's an incoming call.

A CLA2 shaker pad can be plugged into the AmpliCALL20: when the door bell/telephone rings, the shaker pad (optional accessory) will vibrate, making this an ideal solution for heavy sleepers who are expecting important visitors or an important call.

Product Features:

- Extra Loud Ringer (up to 95dB), Adjustable volume
- Different rings for the doorbell and for the telephone ring
- Tone Control (high/medium/low)
- Camera Flash Ring Indicator
- Battery operated 4x AA supplied, Low battery indicator
- 4 alarm options : Off Light + alarm Light + shaker Alarm + shaker
- Wall and Desk Mountable
- Optional power adaptor MP3310 purchase separately

Optional - Geemarc Shaker Ringer Indicator TT-CLA2 RRP AU \$ 50



RRP AU \$199

Diabetic Products



ChillMED MicroCooler -TT-CM-MC

Diabetic Insulin Vial Case

*Safely travel with your diabetes insulin vials, making sure they stay cool for up to 12 hours at a time**

- Stores 2 vials of medication
- 2 mesh elastic pockets on either side for holding other small supplies
- Comes with a re-freezable ice pack
- Up to 12 hours of cooling time*
- Closes securely with zippers
- Product Dimensions: 8" x 2" x 4"

RRP AU \$39



ChillMED Daily Diabetic Bag

Simple design, soft sided daily diabetic case

Offers plenty of storage in an efficiently sized bag

- Two zippered compartments to secure your supplies
- Belt loop on the back of the bag
- Durable nylon construction
- Carrying handle
- Compartments open wide for easy retrieval of supplies
- When empty, the bag measures are 8" x 4.5" x 1"

RRP AU \$39



InsulPak Diabetic Travel Case Bag - TT-IP

Specifically designed with a built-in temperature gauge

Space for your insulin vials, pens, and syringes

- Safe storage of temperature sensitive medicines
- Stays cold between 2° and 10° Celsius for approximately 24-30 hours*
- 24 ounce gel pack included
- Adjustable shoulder strap included
- Window on the exterior of the bag reveals the interior temperature to ensure the safety of your medication
- Insulated centre dividers give you a cold side for medication and non cooled side for supplies
- Ample space to store all supplies
- Can carry up to two weeks of medicines
- Product Dimensions: approx 254mm x 190mm x 76mm
- Main Internal Compartment: 250mm x 170mm x 30mm

RRP AU \$79



Prefilled Syringe Carrying Cases (Pack of 2)


Protect your diabetes insulin syringes

Providing added safety and convenience while travelling

- Includes one White and one Black syringe carrying case
- Works with most insulin syringes (except Monoject)
- Syringe is securely enclosed, Syringes NOT included
- Internal segmented ribbing holds plunger in your correct preset position
- Impact resistant polypropylene

RRP AU \$29

Pill Storage Boxes



7 Day Braille Pill Storage Box – TTPB-7DW

7 Slots per week pill storage box

Load up your medicine for 7 days

- Suitable for storing pills or jewellery
- Each compartment is marked with the day and also marked in Braille for the blind
- Made of high quality materials
- 7 compartments covering 7 days
- Each compartment approx.: 3 x 1.5 x 2.5cm (Total: 13 x 3 x 2.5cm)
- Colour: Solid white

RRP AU \$16




7 Day Braille Pill Storage Box – TTPB-7DG

7 Slots per week pill storage box

Load up your medicine for 7 days

- Suitable for storing pills or jewellery
- Each compartment is marked with the day and also marked in Braille for the blind
- Made of high quality materials
- 7 compartments covering 7 days
- Each compartment approx.: 2.5 x 1.75 x 1.5cm (Total: 12.8 x 2.5 x 1.5cm)
- Colour: Clear green

RRP AU \$16



4x7 Day Clear & Blue Pill Storage Box – TTPB-4x7D

4 Slots per day for 7 days pill storage box

Load up your medicine for 7 days

- Suitable for storing pills or jewellery
- Each compartment is marked with the day
- Made of high quality materials
- 4 compartments per day covering 7 days
- Each compartment approx.: 2.3 x 2.3 x 2cm (Total: 10 x 2.3 x 2cm per day Total box 18 x 11.5 x 2.5)
- Colour: Clear white and blue

RRP AU \$25



MedCenter Traveler XL – TTPB-MCTRAVXL

Large Capacity Weekly Nylon Organizer and Travel Pack

Load up your medicine for 7 days

- VELCRO FLAP allows for quick access to Pill Boxes
- HOLDS UP TO SEVEN LARGE CAPACITY PILL BOXES for up to one full week of medication
- COMPACT SIZE saves space in desk drawer, carry-on luggage or suitcase
- INCLUDES SEVEN LARGE CAPACITY PILL BOXES Numbered with each day of the week
- FOUR LARGE CAPACITY CAVITIES in each included Pill Box

RRP AU \$55

Single Alarm Reminders / Pet Medicine Reminders




7 day, 30 day or 90 day Countdown Timer

Mini alarm clock with set duration for 7 day, 30 day or 90 day countdown depending on model chosen.

- 7 Day, 30 Day or 90 day countdown timer alarm
- Magnetic back will stick on refrigerator or similar materials
- Also has tape which can be used to affix to practically any surface
- Easy to set and use
- When due, alarms sounds 6 times then repeats every 10 seconds until stop pressed or reset
- LCD display
- Volume level approx 30-40db
- 1 x 'LR44' size battery included
- Unit Dimensions: diameter 40 X 17mm height
- 12 Month Limited Warranty by product source (excludes battery)
- Made in Taiwan

TT1-07 7 Day / 1 Week Countdown Timer
TT1-30 30 Day / 1 Month Countdown Timer
TT1-90 90 Day / 3 Month Countdown Timer

RRP AU \$29



Key Ring Timer - TTKR-T

Key ring timer with count-up and count-down function
Small and portable

- Key Ring with count-up and count-down timer
- Count-down from up to 19 hours and 59 minutes
- 2 second pre-warning at 10 and 5 minutes before zero
- Built in Orange LED Torch
- Requires 2 x AAA batteries (not included)
- DOES NOT have storage pill box for storage of medicines
- Breath Alcohol Tester (not guaranteed or warrantied)

RRP AU \$29

TabTimer Products on Display

Prefer to see our products before you buy?

The following organisations have samples of some of our products on display.

Please contact them to see samples of our products before you buy.

Some organisations do not hold stock so in this case please place your order directly with TabTimer.

Please Note: These organisations only have some items available for display and they may require an appointment, please contact them prior to visiting to ensure they have the product you wish to see and to make an appointment. Whilst these organisations may make our products available for display this does not necessarily imply that they endorse or recommend them. It is the responsibility of the individual to make their own assessment as to the suitability of the products.

NSW	
TabTimer Pty Ltd Contact us and we will arrange a convenient Time and place in the Sydney metro area Ph. 1300 TAB TIMER (1300 822 846) Contact: Jonathan Allott	Assistive Technology Australia (AT Aust) [formerly Independent Living Centre] Level 4, Shop 4019 Westpoint Shopping Centre 17 Patrick Street BLACKTOWN NSW 2148 Ph. 1300 452 679 Contact: Ruwani Siriwardene
Independence Australia Health Solutions 6/25-27 Redfern Street WETHERILL PARK NSW 2164 Ph. 1300 788 855 - Contact: Amanda Davison	Parkinson's NSW (for PD NSW Members) Building 17, Macquarie Hospital, 51 Wicks Road, NORTH RYDE NSW 2113 Ph. Info Line; 1800 644 189 – (02) 8051 1900
Mobility 2 You 2 /12 Aston Road ERINA NSW 2250 Ph. (02) 4367 5751	Ansteys Healthcare 70 Bruncker Road BROADMEADOW NSW 2292 Ph. 1300 651 796 - (02) 4040 6161
Mobility Matters Pty Ltd 34 Orient Street BATEMANS BAY NSW 2536 Ph. (02) 4472 1044	Mobility Matters Pty Ltd 117 Carp Street BEGA NSW 2550 Ph. (02) 6492 3888
Epilepsy Action Australia Suite 8, 'Oxford Place', 44-46 Oxford St EPPING NSW 2121 Ph. 1300 37 45 37 - Contact: Gillian Gardner	Calvary Silver Circle - Cessnock 19 Wine Country Dr CESSNOCK NSW 2325 Ph. (02) 4993 9202 - Contact: Jodie Lenton
VIC	
Independence Australia Health Solutions Building 46, Central West Business Pk 9 Ashley Street, WEST FOOTSCRAY VIC 3012 Ph. 1300 788 855	Focus on Lifestyle Factory 3, 17-19 The Concourse COWES VIC 3922 Ph. (03) 5952 3352 - Contact: Fred Lodge
Lyndoch Community Rehabilitation Centre Hopkins Road WARRNAMBOOL VIC 3280 Ph. (03) 5561 9300 - Contact: Kerri Butler	Parkinson's Victoria 587 Canterbury Road Surrey Hills VIC 3127 Ph. 1800 644 189 Or (03) 8809 0400 Contact: Health Information Team

QLD	
LifeTec Queensland (Brisbane) + Mobile Display Level One, Reading Newmarket Cnr Newmarket & Enoggera Roads NEWMARKET QLD 4051 Ph. 1300 885 886 - Contact: Jamie Matveyeff	LifeTec Queensland (Townsville) Domain Central, 103 Duckworth Street GARBUTT QLD 4814 Ph. (07) 4759 5600 - Contact: Judy Jacobsen
Parkinson's QLD U 2, 25 Watland St SPRINGWOOD QLD 4127 Ph. (07) 3209 1588 - Contact: Info Line	Independence Australia Health Solutions 15 Computer Road Yatala QLD 4207 Ph. 1300 788 855 - Contact: Karen Eagle
Epilepsy Queensland Inc Gabba Towers, Level 2 - 411 Vulture St WOOLONGABBA 4102 Ph. 1300 852 853 - Contact: Charlene Mundy	
ACT	
Parkinson's ACT meets monthly @ S.H.O.U.T. Offices Building One, Pearce Community Centre, Collett Place, PEARCE ACT 2607 Ph. (02) 6290 1984	Mobility Matters Pty Ltd 35 Townsville Street FYSHWICK ACT 2609 Ph. 1300 726 906
TAS	
The Independent Living Centre - TAS 275 Wellington Street SOUTH LAUNCESTON TAS 7249 Ph. (03) 6335 9200 - Contact: Pip O'Connor	Independence Australia Health Solutions 225A Kennedy Drive Cambridge, HOBART TAS 7170 Ph. 1300 788 855
SA	
Parkinson's Nurse TAS - Parkinson's Clinic Lower Statton Building, Ground Floor Repatriation Centre, 90 Davey Street, HOBART TAS 7000 Ph.(03) 6166 7763 - Contact: Susan Taylor	The Independent Living Centre - SA + Outreach Disability SA, 11 Blacks Road GILLES PLAINS SA 5086 Ph. 1300 885 886
SA	
Parkinson's SA Inc 23a King William Road UNLEY SA 5061 Ph.(08) 8357 8909 - Contact: Judy White	Independence Australia Health Solutions 3 Duncan Road Dry Creek SA 5094 Ph. 1300 788 855
WA	
Independent Living Centre - WA The Niche, Suite A, 11 Aberdare Road NEDLANDS WA 6009 Ph. 08 6457 7373 - Contact: Lauren Farrell	Parkinson's WA The Niche, Suite B, 11 Aberdare Road NEDLANDS WA 6009 Ph. (08) 6457 7373 - Contact: Angela Sideris
Independence Australia Health Solutions 18-22 Belmont Avenue Rivervale WA 6103 Ph. 1300 788 855	Alzheimer's Aust. WA Assistive Technology Showroom + Mobile 9 Bedbrook Place SHENTON PARK WA 6008 Ph. (08) 9388 2800

About TabTimer

TabTimer™ supplies reminder devices for use by patients on strict medicine regimes, for instance people living with Parkinson's, Alzheimer's, Epilepsy, Diabetes and other conditions that require timely administration of medicines. Some medicines require administration half hourly, hourly, or two hourly and so on. Subsequently, it is often difficult for patients, their relatives, or carers to remember to administer medicines on time.

Reminder devices, which in their simplest form consist of a pill box with a multiple alarm incorporated, aid the process by reminding people and their carers to administer medicines on time. Various styles of reminder devices are available including daily alarms, multiple alarms per day, single or various pill compartments, vibrating alarm watches and clocks. Reminder devices can also be used as a contraceptive pill reminder or pet medication reminder.

The various options will appeal to people depending on the regularity of their medicine schedule.

TabTimer™ helps keep medicines on time™

www.TabTimer.com.au

AU: 1300 TAB TIMER (1300 822 846)

NZ: 0800 TAB TIMER (0800 822 846)

